

# Python Documentation

version

April 15, 2020



# Contents

<b>Welcome to moseq2-pca's documentation!</b>	<b>1</b>
moseq2-pca	1
moseq2_pca package	1
Subpackages	1
moseq2_pca.helpers package	1
Helpers - Data Module	1
Helpers - Wrapper Module	1
moseq2_pca.pca package	2
PCA - Utilities Module	2
moseq2_pca.tests package	5
Subpackages	5
moseq2_pca.tests.integration_tests package	5
Integration Tests - CLI Tests Module	5
Integration Tests - GUI Tests Module	5
moseq2_pca.tests.unit_tests package	5
Unit Tests - PCA Utilities Tests Module	5
Unit Tests - General Utilities Tests Module	5
Unit Tests - Visualization Tests Module	6
CLI Module	6
cli	6
apply-pca	6
clip-scores	7
compute-changepoints	8
train-pca	9
GUI Module	10
Utilities Module	11
Visualization Module	14
<b>Indices and tables</b>	<b>14</b>
<b>Index</b>	<b>15</b>
<b>Python Module Index</b>	<b>21</b>



# Welcome to moseq2-pca's documentation!

## moseq2-pca

### *moseq2\_pca package*

#### *Subpackages*

#### *moseq2\_pca.helpers package*

##### *Helpers - Data Module*

`moseq2_pca.helpers.data.get_pca_yaml_data(pca_yaml)`

Reads PCA yaml file and returns metadata

**Parameters:** **pca\_yaml (str)** (path to pca.yaml)

**Returns:** **use\_fft (bool)** (indicates whether to use FFT) **clean\_params (dict)** (dict of image filtering parameters) **mask\_params (dict)** (dict of mask parameters)) **missing\_data (bool)** (indicates whether to use mask\_params)

`moseq2_pca.helpers.data.load_pcs_for_cp(pca_file_components, config_data)`

Load computed Principal Components.

**Parameters:**

- **pca\_file\_components (str)** (path to pca h5 file to read PCs)
- **config\_data (dict)** (config parameters)

**Returns:** **pca\_components (str)** (path to pca components) **changepoint\_params (dict)** (dict of relevant changepoint parameters) **cluster (dask Cluster)** (Dask Cluster object.) **client (dask Client)** (Dask Client Object) **missing\_data (bool)** (Indicates whether to use mask\_params) **mask\_params (dict)** (Mask parameters to use when computing CPs)

`moseq2_pca.helpers.data.setup_cp_command(input_dir, config_data, output_dir, output_file, output_directory)`

Helper function for changepoints\_wrapper to perform data-path existence checks.

**Parameters:**

- **input\_dir (int)** (path to directory containing all h5+yaml files)
- **config\_data (dict)** (dict of relevant PCA parameters (image filtering etc.))
- **output\_dir (str)** (path to directory to store PCA data)
- **output\_file (str)** (pca model filename)
- **output\_directory (str)** (alternative output\_dir)

**Returns:** **config\_data (dict)** (updated config\_data dict with the proper paths) **pca\_file\_components (str)** (path to trained pca file) **pca\_file\_scores (str)** (path to pca\_scores file) **h5s (list)** (list of relevant pca h5 files) **yamls (list)** (list of relevant pca metadata yaml files) **save\_file (str)** (path to save changepoints)

##### *Helpers - Wrapper Module*

`moseq2_pca.helpers.wrappers.apply_pca_wrapper(input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)`

Wrapper function to obtain PCA Scores.

**Parameters:**

- **input\_dir (int)** (*path to directory containing all h5+yaml files*)
- **config\_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output\_dir (str)** (*path to directory to store PCA data*)
- **output\_file (str)** (*pca model filename*)
- **output\_directory (str)** (*alternative output\_dir*)
- **gui (bool)** (*indicate GUI is running*)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.compute_changepoints_wrapper` (input\_dir, config\_data, output\_dir, output\_file, gui=False, output\_directory=None)

Wrapper function to compute model-free (PCA based) Changepoints.

**Parameters:**

- **input\_dir (int)** (*path to directory containing all h5+yaml files*)
- **config\_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output\_dir (str)** (*path to directory to store PCA data*)
- **output\_file (str)** (*pca model filename*)
- **output\_directory (str)** (*alternative output\_dir*)
- **gui (bool)** (*indicate GUI is running*)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.train_pca_wrapper` (input\_dir, config\_data, output\_dir, output\_file, output\_directory=None, gui=False)

Wrapper function to train PCA.

**Parameters:**

- **input\_dir (int)** (*path to directory containing all h5+yaml files*)
- **config\_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output\_dir (str)** (*path to directory to store PCA data*)
- **output\_file (str)** (*pca model filename*)
- **output\_directory (str)** (*alternative output\_dir*)
- **gui (bool)** (*indicate GUI is running*)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

## ***moseq2\_pca.pca package***

### ***PCA - Utilities Module***

`moseq2_pca.pca.util.apply_pca_dask` (pca\_components, h5s, yamls, use\_fft, clean\_params, save\_file, chunk\_size, mask\_params, missing\_data, client, fps=30, gui=False)

“Apply” trained PCA on input frame data to obtain PCA Scores using Distributed Dask cluster.

**Parameters:**

- **pca\_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use\_fft (bool)** (indicate whether to use 2D-FFT)
- **clean\_params (dict)** (dictionary containing filtering options)
- **save\_file (str)** (path to pca\_scores filename to save)
- **chunk\_size (int)** (size of chunks to process)
- **mask\_params (dict)** (dictionary of masking parameters (if missing data))
- **missing\_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.apply_pca_local` (pca\_components, h5s, yamls, use\_fft, clean\_params, save\_file, chunk\_size, mask\_params, missing\_data, fps=30)

“Apply” trained PCA on input frame data to obtain PCA Scores using local cluster/platform.

**Parameters:**

- **pca\_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use\_fft (bool)** (indicate whether to use 2D-FFT)
- **clean\_params (dict)** (dictionary containing filtering options)
- **save\_file (str)** (path to pca\_scores filename to save)
- **chunk\_size (int)** (size of chunks to process)
- **mask\_params (dict)** (dictionary of masking parameters (if missing data))
- **missing\_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.get_changepoints_dask` (changepoint\_params, pca\_components, h5s, yamls, save\_file, chunk\_size, mask\_params, missing\_data, client, fps=30, pca\_scores=None, progress\_bar=False, gui=False)

Computes model-free changepoints using PCs and PC Scores on distributed dask cluster.

**Parameters:**

- **changepoint\_params (dict)** (*dict of changepoint parameters*)
- **pca\_components (np.array)** (*computed principal components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **save\_file (str)** (*path to save changepoint files*)
- **chunk\_size (int)** (*size of chunks to process in dask.*)
- **mask\_params (dict)** (*dict of missing\_data mask parameters.*)
- **missing\_data (bool)** (*indicate whether to use mask\_params*)
- **client (dask Client)** (*initialized Dask Client object*)
- **fps (int)** (*frames per second*)
- **pca\_scores (np.array)** (*computed principal component scores*)
- **progress\_bar (bool)** (*display progress bar*)
- **gui (bool)** (*indicate GUI use*)

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.mask_data` (*original\_data, mask, new\_data*)  
Create a mask subregion given a boolean mask if missing data flag is used.

**Parameters:**

- **original\_data (3d np.ndarray)** (*input frames*)
- **mask (3d boolean np.ndarray)** (*mask array*)
- **new\_data (3d np.ndarray)** (*frames to use*)

**Returns:** output (3d np.ndarray)

**Return type:** masked data array

`moseq2_pca.pca.util.train_pca_dask` (*dask\_array, clean\_params, use\_fft, rank, cluster\_type, client, workers, cache, mask=None, iters=10, recon\_pcs=10, min\_height=10, max\_height=100*)  
Train PCA using dask arrays.

**Parameters:**

- **dask\_array (dask array)** (*chunked frames to train PCA*)
- **clean\_params (dict)** (*dictionary containing filtering parameters*)
- **use\_fft (bool)** (*indicates whether to use 2d-FFT on images.*)
- **rank (int)** (*Matrix rank to use*)
- **cluster\_type (str)** (*indicates which cluster to use.*)
- **client (Dask.Client)** (*client object to execute dask operations*)
- **workers (int)** (*number of dask workers*)
- **cache (str)** (*path to cache directory*)
- **mask (dask array)** (*dask array of masked data if missing\_data parameter==True*)
- **iters (int)** (*number of SVD iterations*)
- **recon\_pcs (int)** (*number of PCs to reconstruct. (if missing\_data = True)*)
- **min\_height (int)** (*minimum mouse height from floor in (mm)*)
- **max\_height (int)** (*maximum mouse height from floor in (mm)*)

**Returns:** output\_dict (dict)

**Return type:** dictionary containing PCA training results.



## ***moseq2\_pca.tests package***

### ***Subpackages***

## ***moseq2\_pca.tests.integration\_tests package***

### ***Integration Tests - CLI Tests Module***

```
class moseq2_pca.tests.integration_tests.test_cli.TestCli (methodName='runTest')
    Bases: unittest.case.TestCase

    test_apply_pca ()

    test_clip_scores ()

    test_compute_changepoints ()

    test_train_pca ()
```

### ***Integration Tests - GUI Tests Module***

```
class moseq2_pca.tests.integration_tests.test_gui.TestGUI (methodName='runTest')
    Bases: unittest.case.TestCase

    test_apply_pca_command ()

    test_compute_changepoints_command ()

    test_train_pca_command ()
```

## ***moseq2\_pca.tests.unit\_tests package***

### ***Unit Tests - PCA Utilities Tests Module***

```
class moseq2_pca.tests.unit_tests.test_pca_util.TestPCAUtils (methodName='runTest')
    Bases: unittest.case.TestCase

    test_mask_data ()

    test_train_pca_dask ()
```

### ***Unit Tests - General Utilities Tests Module***

```
class moseq2_pca.tests.unit_tests.test_util.TestUtils (methodName='runTest')
    Bases: unittest.case.TestCase

    test_clean_frames ()

    test_gauss_smooth ()

    test_gaussian_kernel1d ()

    test_get_changepoints ()

    test_get_metadata_path ()
```

Welcome to moseq2-pca's documentation!

```
test_get_rps ()
test_get_rsp_dask ()
test_get_timestamp_path ()
test_initialize_dask ()
test_insert_nans ()
test_read_yaml ()
test_recursive_find_h5s ()
test_select_strel ()
```

### Unit Tests - Visualization Tests Module

```
class moseq2_pca.tests.unit_tests.test_viz.TestViz (methodName='runTest')
    Bases: unittest.case.TestCase

    changepoint_dist ()

    test_display_components ()
        cmap = 'gray' im_size = int(np.sqrt(components.shape[1])) plotv = components.reshape((-1, im_size, im_size))
        plotv = skimage.util.montage(plotv)
        plt.switch_backend('agg')
        fig, ax = plt.subplots(1, 1, figsize=(10, 10)) plt.imshow(plotv, cmap=cmap) plt.xticks([]) plt.yticks([])

    test_scee_plot ()
```

## CLI Module

### cli

```
cli [OPTIONS] COMMAND [ARGS]...
```

### apply-pca

```
cli apply-pca [OPTIONS]
```

#### Options

```
-i, --input-dir <input_dir>
    Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

--cluster-type <cluster_type>
    Cluster type [default: local]

    Options:  local|slurm|nodask

-o, --output-dir <output_dir>
    Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca]

--output-file <output_file>
    Name of h5 file for storing pca results [default: pca_scores]

--h5-path <h5_path>
    Path to data in h5 files [default: /frames]

--h5-mask-path <h5_mask_path>
```

Path to log-likelihood mask in h5 files [default: /frames\_mask]

**--pca-path** <pca\_path>  
Path to pca components [default: /components]

**--pca-file** <pca\_file>  
Path to PCA results

**--chunk-size** <chunk\_size>  
Number of frames per chunk [default: 4000]

**--fill-gaps** <fill\_gaps>  
Fill dropped frames with nans [default: True]

**--fps** <fps>  
Fps (only used if no timestamps found) [default: 30]

**--detrend-window** <detrend\_window>  
Length of detrend window (in seconds, 0 for no detrending) [default: 0]

**--config-file** <config\_file>  
Path to configuration file

**-d, --dask-cache-path** <dask\_cache\_path>  
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2\_pca]

**-q, --queue** <queue>  
Cluster queue/partition for submitting jobs [default: debug]

**-n, --nworkers** <nworkers>  
Number of workers [default: 10]

**-c, --cores** <cores>  
Number of cores per worker [default: 1]

**-p, --processes** <processes>  
Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>  
RAM usage per workers [default: 15GB]

**-w, --wall-time** <wall\_time>  
Wall time for workers [default: 06:00:00]

**--timeout** <timeout>  
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

## clip-scores

Clips PCA scores from the beginning or end

### Args:

**pca\_file** (string): Path to PCA scores  
**clip\_samples** (int): number of samples to clip from beginning or end  
**from\_end** (bool): if true clip from end rather than beginning

Note that scores are modified *in place*.

```
cli clip-scores [OPTIONS] PCA_FILE CLIP_SAMPLES
```

### Options

**--from-end**  
[default: False]

### Arguments

**PCA\_FILE**  
Required argument

**CLIP\_SAMPLES**  
Required argument

## **compute-changepoints**

```
cli compute-changepoints [OPTIONS]
```

### Options

**-i, --input-dir** <input\_dir>  
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

**-o, --output-dir** <output\_dir>  
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/\_pca/]

**--output-file** <output\_file>  
Name of h5 file for storing pca results [default: changepoints]

**--cluster-type** <cluster\_type>  
Cluster type [default: local]

**Options:** local|slurm

**--pca-file-components** <pca\_file\_components>  
Path to PCA components

**--pca-file-scores** <pca\_file\_scores>  
Path to PCA results

**--pca-path** <pca\_path>  
Path to pca components [default: /components]

**--neighbors** <neighbors>  
Neighbors to use for peak identification [default: 1]

**--threshold** <threshold>  
Peak threshold to use for changepoints [default: 0.5]

**-k, --klags** <klags>  
Lag to use for derivative calculation [default: 6]

**-s, --sigma** <sigma>  
Standard deviation of gaussian smoothing filter [default: 3.5]

**-d, --dims** <dims>  
Number of random projections to use [default: 300]

**--fps** <fps>  
Fps (only used if no timestamps found) [default: 30]

**--h5-path** <h5\_path>  
Path to data in h5 files [default: /frames]

**--h5-mask-path** <h5\_mask\_path>  
Path to log-likelihood mask in h5 files [default: /frames\_mask]

**--chunk-size** <chunk\_size>  
Number of frames per chunk [default: 4000]

**--config-file** <config\_file>  
Path to configuration file

**--dask-cache-path** <dask\_cache\_path>  
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2\_pca]

**--visualize-results** <visualize\_results>  
Visualize results [default: True]

**-q, --queue** <queue>  
Cluster queue/partition for submitting jobs [default: debug]

**-n, --nworkers** <nworkers>  
Number of workers [default: 10]

**-c, --cores** <cores>

Welcome to moseq2-pca's documentation!

Number of cores per worker [default: 1]

**-p, --processes** <processes>

Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>

RAM usage per workers [default: 15GB]

**-w, --wall-time** <wall\_time>

Wall time for workers [default: 06:00:00]

**--timeout** <timeout>

Time to wait for workers to initialize before proceeding (minutes) [default: 5]

## ***train-pca***

```
cli train-pca [OPTIONS]
```

### Options

**-i, --input-dir** <input\_dir>

Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

**--cluster-type** <cluster\_type>

Cluster type [default: local]

**Options:** local|slurm

**-o, --output-dir** <output\_dir>

Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/\_pca]

**--gaussfilter-space** <gaussfilter\_space>

Spatial filter for data (Gaussian) [default: 1.5, 1]

**--gaussfilter-time** <gaussfilter\_time>

Temporal filter for data (Gaussian) [default: 0]

**--medfilter-space** <medfilter\_space>

Median spatial filter [default: 0]

**--medfilter-time** <medfilter\_time>

Median temporal filter [default: 0]

**--missing-data**

Use missing data PCA [default: False]

**--missing-data-iters** <missing\_data\_iters>

Missing data PCA iterations [default: 10]

**--mask-threshold** <mask\_threshold>

Threshold for mask (missing data only) [default: -16]

**--mask-height-threshold** <mask\_height\_threshold>

Threshold for mask based on floor height [default: 5]

**--min-height** <min\_height>

Min mouse height from floor (mm) [default: 10]

**--max-height** <max\_height>

Max mouse height from floor (mm) [default: 100]

**--tailfilter-size** <tailfilter\_size>

Tail filter size [default: 9, 9]

**--tailfilter-shape** <tailfilter\_shape>

Tail filter shape [default: ellipse]

**--use-fft**

Use 2D fft [default: False]

**--recon-pcs** <recon\_pcs>

Welcome to moseq2-pca's documentation!

Number of PCs to use for missing data reconstruction [default: 10]

**--rank** <rank>  
Rank for compressed SVD (generally>>nPCS) [default: 50]

**--output-file** <output\_file>  
Name of h5 file for storing pca results [default: pca]

**--chunk-size** <chunk\_size>  
Number of frames per chunk [default: 4000]

**--visualize-results** <visualize\_results>  
Visualize results [default: True]

**--config-file** <config\_file>  
Path to configuration file

**-d, --dask-cache-path** <dask\_cache\_path>  
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2\_pca]

**--local-processes** <local\_processes>  
Use processes with local scheduler [default: True]

**-q, --queue** <queue>  
Cluster queue/partition for submitting jobs [default: debug]

**-n, --nworkers** <nworkers>  
Number of workers [default: 10]

**-c, --cores** <cores>  
Number of cores per worker [default: 1]

**-p, --processes** <processes>  
Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>  
Total RAM usage per worker [default: 15GB]

**-w, --wall-time** <wall\_time>  
Wall time for workers [default: 06:00:00]

**--timeout** <timeout>  
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

## GUI Module

`moseq2_pca.gui.apply_pca_command` (input\_dir, index\_file, config\_file, output\_dir, output\_file, output\_directory=None)

Compute PCA Scores given trained PCA using Jupyter Notebook.

**Parameters:**

- **input\_dir (str)** (path to directory containing training data)
- **index\_file (str)** (path to index file.)
- **config\_file (str)** (path to config file)
- **output\_dir (str)** (path to output pca directory)
- **output\_file (str)** (name of output pca file.)
- **output\_directory (str)** (alternative output directory path)

**Returns:** (str)

**Return type:** success string.

`moseq2_pca.gui.compute_changepoints_command` (input\_dir, config\_file, output\_dir, output\_file, output\_directory=None)

Compute Changepoint distribution using Jupyter Notebook.

**Parameters:**

- **input\_dir (str)** (*path to directory containing training data*)
- **config\_file (str)** (*path to config file*)
- **output\_dir (str)** (*path to output pca directory*)
- **output\_file (str)** (*name of output pca file.*)
- **output\_directory (str)** (*alternative output directory path*)

**Returns:** (str)

**Return type:** success string.

`moseq2_pca.gui.train_pca_command` (input\_dir, config\_file, output\_dir, output\_file, output\_directory=None)  
Train PCA through Jupyter notebook, and updates config file.

**Parameters:**

- **input\_dir (str)** (*path to directory containing training data*)
- **config\_file (str)** (*path to config file*)
- **output\_dir (str)** (*path to output pca directory*)
- **output\_file (str)** (*name of output pca file.*)
- **output\_directory (str)** (*alternative output directory path*)

**Returns:**

**Return type:** None

## Utilities Module

`moseq2_pca.util.clean_frames` (frames, medfilter\_space=None, gaussfilter\_space=None, medfilter\_time=None, gaussfilter\_time=None, detrend\_time=None, tailfilter=None, tail\_threshold=5)  
Filters spatial/temporal noise from frames using Median and Gaussian filters, given kernel sizes for each respective requested filter.

**Parameters:**

- **frames (3D numpy array)** (*frames to filter.*)
- **medfilter\_space (list)** (*median spatial filter kernel.*)
- **gaussfilter\_space (list)** (*gaussian spatial filter kernel.*)
- **medfilter\_time (list)** (*median temporal filter.*)
- **gaussfilter\_time (list)** (*gaussian temporal filter.*)
- **detrend\_time (int)** (*number of frames to lag for.*)
- **tailfilter (int)** (*size of tail-filter kernel.*)
- **tail\_threshold (int)** (*threshold value to use for tail filtering*)

**Returns:** out (3D numpy array)

**Return type:** filtered frames.

`moseq2_pca.util.command_with_config` (config\_file\_param\_name)

`moseq2_pca.util.gauss_smooth` (signal, win\_length=None, sig=1.5, kernel=None)  
Perform Gaussian Smoothing on a 1D signal.

**Parameters:**

- **signal (1d numpy array)** (*signal to perform smoothing*)
- **win\_length (int)** (*window\_size for gaussian kernel filter*)
- **sig (float)** (*variance of 1d gaussian kernel.*)
- **kernel (tuple)** (*kernel size to use for smoothing*)

**Returns:** result (1d numpy array)

**Return type:** smoothed signal

Welcome to moseq2-pca's documentation!

`moseq2_pca.util.gaussian_kernel1d` (`n=None`, `sig=3`)

Get 1D gaussian kernel.

**Parameters:**

- **n (int)** (*number of points to use.*)
- **sig (int)** (*variance of kernel to use.*)

**Returns:** **kernel (1d array)**

**Return type:** 1D numpy kernel.

`moseq2_pca.util.get_changepoints` (`scores`, `k=5`, `sigma=3`, `peak_height=0.5`, `peak_neighbors=1`, `baseline=True`, `timestamps=None`)

Compute changepoints distribution and CP Curve.

**Parameters:**

- **scores (3D numpy array)** (*nframes \* r \* c*)
- **k (int)** (*klags - Lag to use for derivative calculation.*)
- **sigma (int)** (*Standard deviation of gaussian smoothing filter.*)
- **peak\_height (float)** (*user-defined peak Changepoint length.*)
- **peak\_neighbors (int)** (*number of peaks in the CP curve.*)
- **baseline (bool)** (*normalize data.*)
- **timestamps (array)** (*loaded timestamps.*)

**Returns:** **cps (numpy array)** (*array of values for CP curve*) **normed\_df (numpy array)** (*array of values for bar plot*)

`moseq2_pca.util.get_metadata_path` (`h5file`)

Return path within h5 file that contains the kinect extraction metadata.

**Parameters:** **h5file (str)** (*path to h5 file.*)

**Returns:** **(str)**

**Return type:** path to acquisition metadata within h5 file.

`moseq2_pca.util.get_rps` (`frames`, `rps=600`, `normalize=True`)

Get random projections of frames.

**Parameters:**

- **frames (2D or 3D numpy array)** (*Frames to get dimensions from.*)
- **rps (int)** (*Number of random projections.*)
- **normalize (bool)** (*indicates whether to normalize frames.*)

**Returns:** **rproj (2D or 3D numpy array)**

**Return type:** Computed random projections with same shape as frames

`moseq2_pca.util.get_timestamp_path` (`h5file`)

Return path within h5 file that contains the kinect timestamps

**Parameters:** **h5file (str)** (*path to h5 file.*)

**Returns:** **(str)**

**Return type:** path to metadata timestamps within h5 file

`moseq2_pca.util.initialize_dask` (`nworkers=50`, `processes=1`, `memory='4GB'`, `cores=1`, `wall_time='01:00:00'`, `queue='debug'`, `local_processes=False`, `cluster_type='local'`, `scheduler='distributed'`, `timeout=10`, `cache_path='/Users/aymanzeine/moseq2_pca'`, `**kwargs`)

Initialize dask client, cluster, workers, etc.



**Parameters:**

- **nworkers (int)** (*number of dask workers to initialize*)
- **processes (int)** (*number of processes per worker*)
- **memory (str)** (*amount of memory to allocate to dask cluster*)
- **cores (int)** (*number of cores to use.*)
- **wall\_time (str)** (*amount of time to allow program to run*)
- **queue (str)** (*logging mode*)
- **local\_processes (bool)** (*indicate whether the processes are local*)
- **cluster\_type (str)** (*indicate what cluster to use*)
- **scheduler (str)** (*indicate what scheduler to use*)
- **timeout (int)** (*number of worker timeouts to allow*)
- **cache\_path (str or Pathlike)** (*path to store cached data*)
- **kwargs** (*extra keyword arguments*)

**Returns:** **client (dask Client)** (*initialized Client*) **cluster (dask Cluster)** (*initialized Cluster*) **workers (dask Workers)** (*intialized workers*) **cache (dask Chest)** (*initialized Chest (cache) object*)

`moseq2_pca.util.insert_nans` (timestamps, data, fps=30)  
Fills NaN values with 0 in timestamps.

**Parameters:**

- **timestamps (1D array)** (*timestamp time-strs*)
- **data (1D array)** (*timestamp values*)
- **fps (int)** (*frames per second*)

**Returns:** **filled\_data (1D array)** (*filled missing timestamp values.*) **data\_idx (1D array)** (*indices of inserted 0s*) **filled\_timestamps (1D array)** (*filled timestamp-strs*)

`moseq2_pca.util.read_yaml` (yaml\_file)  
Reads yaml file and returns dictionary representation of file contents.

**Parameters:** **yaml\_file (str)** (*path to yaml file*)

**Returns:** **return\_dict (dict)**

**Return type:** dict of yaml file contents

`moseq2_pca.util.recursive_find_h5s` (root\_dir='/Users/aymanzeine/Desktop/moseq/moseq2-pca/docs', ext='.h5', yaml\_string='{}.yaml')  
Recursively find h5 files, along with yaml files with the same basename

**Parameters:**

- **root\_dir (str or os.Pathlike)** (*path to directory to start recursive search*)
- **ext (str)** (*extension to search for, e.g. .h5*)
- **yaml\_string (str)** (*a format to use to name yaml files*)

**Returns:** **h5s (list)** (*list of h5 file paths*) **dicts (list)** (*list of metadata file paths*) **yamls (list)** (*list of yaml file paths*)

`moseq2_pca.util.recursively_load_dict_contents_from_group` (h5file, path)  
Reads all contents from h5 and returns them in a nested dict object.

**Parameters:**

- **h5file (str)** (*path to h5 file*)
- **path (str)** (*path to group within h5 file*)

**Returns:** **ans (dict)**

**Return type:** dictionary of all h5 group contents

`moseq2_pca.util.select_strel` (string='e', size=10, 10)  
Selects Structuring Element Shape

**Parameters:**

- **string (str)** (*e for Ellipse, r for Rectangle*)
- **size (tuple)** (*size of StructuringElement*)

**Returns:** **strel (cv2.StructuringElement)****Return type:** returned StructuringElement with specified size.`moseq2_pca.util.shutdown_dask (scheduler)`

Graceful

shutdown

dask

scheduler.

source:

<https://github.com/dask/distributed/issues/1703#issuecomment-361291492>**Parameters:** **scheduler (dask Scheduler)** (*scheduler to shutdown.*)**Returns:****Return type:** None

## Visualization Module

`moseq2_pca.viz.changepoint_dist (cps, headless=False)`

Creates bar plot describing computed Changepoint Distribution.

**Parameters:**

- **cps (np.ndarray)** (*changepoints to graph*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)`moseq2_pca.viz.display_components (components, cmap='gray', headless=False)`

Creates grid of computed Principal Components.

**Parameters:**

- **components (np.ndarray)** (*components to graph*)
- **cmap (str)** (*color map to use*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)`moseq2_pca.viz.scrree_plot (explained_variance_ratio, headless=False)`

Creates Scree plot describing principal components.

**Parameters:**

- **explained\_variance\_ratio (np.array)** (*explained variance ratio of each principal component*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)****Return type:** figure to save/graph

## Indices and tables

- `genindex`
- `modindex`
- `search`

# Index

## Symbols

			--gaussfilter-time <gaussfilter_time>	cli-train-pca command line option
			--h5-mask-path <h5_mask_path>	cli-apply-pca line option command
--chunk-size <chunk_size>	cli-apply-pca line option	command		cli-compute-changepoints command line option
	cli-compute-changepoints command line option		--h5-path <h5_path>	cli-apply-pca line option command
	cli-train-pca line option	command		cli-compute-changepoints command line option
--cluster-type <cluster_type>	cli-apply-pca line option	command	--input-dir <input_dir>	cli-apply-pca line option command
	cli-compute-changepoints command line option			cli-compute-changepoints command line option
	cli-train-pca line option	command		cli-train-pca line option command
--config-file <config_file>	cli-apply-pca line option	command	--klags <klags>	cli-compute-changepoints command line option
	cli-compute-changepoints command line option		local-processes <local_processes>	cli-train-pca command line option
	cli-train-pca line option	command	--mask-height-threshold <mask_height_threshold>	cli-train-pca command line option
--cores <cores>	cli-apply-pca line option	command		
	cli-compute-changepoints command line option		--mask-threshold <mask_threshold>	cli-train-pca command line option
	cli-train-pca line option	command	--max-height <max_height>	cli-train-pca command line option
--dask-cache-path <dask_cache_path>	cli-apply-pca line option	command	--medfilter-space <medfilter_space>	cli-train-pca command line option
	cli-compute-changepoints command line option		--medfilter-time <medfilter_time>	cli-train-pca command line option
	cli-train-pca line option	command		
--detrend-window <detrend_window>	cli-apply-pca command line option		--memory <memory>	cli-apply-pca line option command
	cli-compute-changepoints command line option			cli-compute-changepoints command line option
--dims <dims>	cli-apply-pca command line option			cli-train-pca line option command
--fill-gaps <fill_gaps>	cli-apply-pca command line option		--min-height <min_height>	cli-train-pca command line option
--fps <fps>	cli-apply-pca line option	command	--missing-data	cli-train-pca line option command
	cli-compute-changepoints command line option		--missing-data-itsers <missing_data_itsers>	cli-train-pca command line option
--from-end	cli-clip-scores line option	command		
--gaussfilter-space <gaussfilter_space>	cli-train-pca command line option		--neighbors <neighbors>	cli-compute-changepoints command line option

--nworkers <nworkers>	cli-apply-pca line option	command	--threshold <threshold>	cli-compute-change points command line option
	cli-compute-change points command line option		--timeout <timeout>	cli-apply-pca line option
	cli-train-pca line option	command		cli-compute-change points command line option
--output-dir <output_dir>	cli-apply-pca line option	command		cli-train-pca line option
	cli-compute-change points command line option		--use-fft	cli-train-pca line option
	cli-train-pca line option	command	--visualize-results <visualize_results>	cli-compute-change points command line option
--output-file <output_file>	cli-apply-pca line option	command		cli-train-pca line option
	cli-compute-change points command line option		--wall-time <wall_time>	cli-apply-pca line option
	cli-train-pca line option	command		cli-compute-change points command line option
--pca-file <pca_file>	cli-apply-pca command line option			cli-train-pca line option
-file-components <pca_file_components>	cli-compute-change points command line option			cli-apply-pca line option
-pca-file-scores <pca_file_scores>	cli-compute-change points command line option			cli-compute-change points command line option
--pca-path <pca_path>	cli-apply-pca line option	command		cli-train-pca line option
	cli-compute-change points command line option		-d	cli-apply-pca line option
--processes <processes>	cli-apply-pca line option	command		cli-compute-change points command line option
	cli-compute-change points command line option			cli-train-pca line option
	cli-train-pca line option	command	-i	cli-apply-pca line option
--queue <queue>	cli-apply-pca line option	command		cli-compute-change points command line option
	cli-compute-change points command line option			cli-train-pca line option
	cli-train-pca line option	command	-k	cli-compute-change points command line option
--rank <rank>	cli-train-pca line option	command	-m	cli-apply-pca line option
--recon-pcs <recon_pcs>	cli-train-pca command line option			cli-compute-change points command line option
--sigma <sigma>	cli-compute-change points command line option			cli-train-pca line option
--tailfilter-shape <tailfilter_shape>	cli-train-pca command line option		-n	cli-apply-pca line option
	cli-train-pca command line option			cli-compute-change points command line option
--tailfilter-size <tailfilter_size>	cli-train-pca command line option			cli-train-pca line option

-o	cli-apply-pca	command line option	--h5-mask-path <h5_mask_path>
	cli-compute-changepoints	command line option	--h5-path <h5_path>
	cli-train-pca	command line option	--input-dir <input_dir>
-p	cli-apply-pca	command line option	--memory <memory>
	cli-compute-changepoints	command line option	--nworkers <nworkers>
	cli-train-pca	command line option	--output-dir <output_dir>
-q	cli-apply-pca	command line option	--output-file <output_file>
	cli-compute-changepoints	command line option	--pca-file <pca_file>
	cli-train-pca	command line option	--pca-path <pca_path>
-s	cli-apply-pca	command line option	--processes <processes>
	cli-compute-changepoints	command line option	--queue <queue>
	cli-train-pca	command line option	--timeout <timeout>
-w	cli-apply-pca	command line option	--wall-time <wall_time>
	cli-compute-changepoints	command line option	-c
	cli-train-pca	command line option	-d
			-i
			-m
			-n
			-o
			-p
			-q
			-w

## A

apply\_pca\_command() (in module moseq2\_pca.gui)

apply\_pca\_dask() (in module moseq2\_pca.pca.util)

apply\_pca\_local() (in module moseq2\_pca.pca.util)

apply\_pca\_wrapper() (in module moseq2\_pca.helpers.wrappers)

## C

changepoint\_dist() (in module moseq2\_pca.viz)  
(moseq2\_pca.tests.unit\_tests.test\_viz.TestViz method)

clean\_frames() (in module moseq2\_pca.util)

### cli-apply-pca command line option

--chunk-size <chunk\_size>  
--cluster-type <cluster\_type>  
--config-file <config\_file>  
--cores <cores>  
--dask-cache-path <dask\_cache\_path>  
--detrend-window <detrend\_window>  
--fill-gaps <fill\_gaps>  
--fps <fps>

### cli-clip-scores command line option

--from-end

CLIP\_SAMPLES

PCA\_FILE

### cli-compute-changepoints command line option

--chunk-size <chunk\_size>  
--cluster-type <cluster\_type>  
--config-file <config\_file>  
--cores <cores>  
--dask-cache-path <dask\_cache\_path>  
--dims <dims>  
--fps <fps>  
--h5-mask-path <h5\_mask\_path>  
--h5-path <h5\_path>  
--input-dir <input\_dir>  
--klags <klags>  
--memory <memory>  
--neighbors <neighbors>  
--nworkers <nworkers>  
--output-dir <output\_dir>

```

--output-file <output_file>
--pca-file-components <pca_file_components>
--pca-file-scores <pca_file_scores>
--pca-path <pca_path>
--processes <processes>
--queue <queue>
--sigma <sigma>
--threshold <threshold>
--timeout <timeout>
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-k
-m
-n
-o
-p
-q
-s
-w

```

#### cli-train-pca command line option

```

--chunk-size <chunk_size>
--cluster-type <cluster_type>
--config-file <config_file>
--cores <cores>
--dask-cache-path <dask_cache_path>
--gaussfilter-space <gaussfilter_space>
--gaussfilter-time <gaussfilter_time>
--input-dir <input_dir>
--local-processes <local_processes>
--mask-height-threshold <mask_height_threshold>
--mask-threshold <mask_threshold>
--max-height <max_height>
--medfilter-space <medfilter_space>
--medfilter-time <medfilter_time>
--memory <memory>
--min-height <min_height>
--missing-data
--missing-data-iters <missing_data_iters>
--nworkers <nworkers>

```

```

--output-dir <output_dir>
--output-file <output_file>
--processes <processes>
--queue <queue>
--rank <rank>
--recon-pcs <recon_pcs>
--tailfilter-shape <tailfilter_shape>
--tailfilter-size <tailfilter_size>
--timeout <timeout>
--use-fft
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-m
-n
-o
-p
-q
-w

```

#### CLIP\_SAMPLES

cli-clip-scores command line option

```

command_with_config() (in module moseq2_pca.util)
compute_changepoints_command() (in module
moseq2_pca.gui)
compute_changepoints_wrapper() (in module
moseq2_pca.helpers.wrappers)

```

#### D

display\_components() (in module moseq2\_pca.viz)

#### G

```

gauss_smooth() (in module moseq2_pca.util)
gaussian_kernel1d() (in module moseq2_pca.util)
get_changepoints() (in module moseq2_pca.util)
get_changepoints_dask() (in module
moseq2_pca.pca.util)
get_metadata_path() (in module moseq2_pca.util)
get_pca_yaml_data() (in module
moseq2_pca.helpers.data)
get_rps() (in module moseq2_pca.util)
get_timestamp_path() (in module moseq2_pca.util)

```

## I

[initialize\\_dask\(\)](#) (in module [moseq2\\_pca.util](#))  
[insert\\_nans\(\)](#) (in module [moseq2\\_pca.util](#))

## L

[load\\_pcs\\_for\\_cp\(\)](#) (in module [moseq2\\_pca.helpers.data](#))

## M

[mask\\_data\(\)](#) (in module [moseq2\\_pca.pca.util](#))

### module

[moseq2\\_pca.gui](#)  
[moseq2\\_pca.helpers.data](#)  
[moseq2\\_pca.helpers.wrappers](#)  
[moseq2\\_pca.pca.util](#)  
[moseq2\\_pca.tests.integration\\_tests.test\\_cli](#)  
[moseq2\\_pca.tests.integration\\_tests.test\\_gui](#)  
[moseq2\\_pca.tests.unit\\_tests.test\\_pca\\_util](#)  
[moseq2\\_pca.tests.unit\\_tests.test\\_util](#)  
[moseq2\\_pca.tests.unit\\_tests.test\\_viz](#)  
[moseq2\\_pca.util](#)  
[moseq2\\_pca.viz](#)

### [moseq2\\_pca.gui](#)

module

### [moseq2\\_pca.helpers.data](#)

module

### [moseq2\\_pca.helpers.wrappers](#)

module

### [moseq2\\_pca.pca.util](#)

module

### [moseq2\\_pca.tests.integration\\_tests.test\\_cli](#)

module

### [moseq2\\_pca.tests.integration\\_tests.test\\_gui](#)

module

### [moseq2\\_pca.tests.unit\\_tests.test\\_pca\\_util](#)

module

### [moseq2\\_pca.tests.unit\\_tests.test\\_util](#)

module

### [moseq2\\_pca.tests.unit\\_tests.test\\_viz](#)

module

### [moseq2\\_pca.util](#)

module

### [moseq2\\_pca.viz](#)

module

## P

### PCA\_FILE

[cli-clip-scores](#) command line option

## R

[read\\_yaml\(\)](#) (in module [moseq2\\_pca.util](#))  
[recursive\\_find\\_h5s\(\)](#) (in module [moseq2\\_pca.util](#))  
[recursively\\_load\\_dict\\_contents\\_from\\_group\(\)](#) (in module [moseq2\\_pca.util](#))

## S

[screed\\_plot\(\)](#) (in module [moseq2\\_pca.viz](#))  
[select\\_strel\(\)](#) (in module [moseq2\\_pca.util](#))  
[setup\\_cp\\_command\(\)](#) (in module [moseq2\\_pca.helpers.data](#))  
[shutdown\\_dask\(\)](#) (in module [moseq2\\_pca.util](#))

## T

[test\\_apply\\_pca\(\)](#)  
([moseq2\\_pca.tests.integration\\_tests.test\\_cli.TestCli](#) method)  
[test\\_apply\\_pca\\_command\(\)](#)  
([moseq2\\_pca.tests.integration\\_tests.test\\_gui.TestGUI](#) method)  
[test\\_clean\\_frames\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_util.TestUtils](#) method)  
[test\\_clip\\_scores\(\)](#)  
([moseq2\\_pca.tests.integration\\_tests.test\\_cli.TestCli](#) method)  
[test\\_compute\\_changepoints\(\)](#)  
([moseq2\\_pca.tests.integration\\_tests.test\\_cli.TestCli](#) method)  
[test\\_compute\\_changepoints\\_command\(\)](#)  
([moseq2\\_pca.tests.integration\\_tests.test\\_gui.TestGUI](#) method)  
[test\\_display\\_components\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_viz.TestViz](#) method)  
[test\\_gauss\\_smooth\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_util.TestUtils](#) method)  
[test\\_gaussian\\_kernel1d\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_util.TestUtils](#) method)  
[test\\_get\\_changepoints\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_util.TestUtils](#) method)  
[test\\_get\\_metadata\\_path\(\)](#)  
([moseq2\\_pca.tests.unit\\_tests.test\\_util.TestUtils](#) method)

[test\\_get\\_rps\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_get\\_rsp\\_dask\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_get\\_timestamp\\_path\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_initialize\\_dask\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_insert\\_nans\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_mask\\_data\(\)](#) (moseq2\_pca.tests.unit\_tests.test\_pca\_util.TestPCAUtils method)

[test\\_read\\_yaml\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_recursive\\_find\\_h5s\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_scree\\_plot\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_viz.TestViz method)

[test\\_select\\_strel\(\)](#)  
 (moseq2\_pca.tests.unit\_tests.test\_util.TestUtils method)

[test\\_train\\_pca\(\)](#)  
 (moseq2\_pca.tests.integration\_tests.test\_cli.TestCli method)

[test\\_train\\_pca\\_command\(\)](#)  
 (moseq2\_pca.tests.integration\_tests.test\_gui.TestGUI method)

[test\\_train\\_pca\\_dask\(\)](#) (moseq2\_pca.tests.unit\_tests.test\_pca\_util.TestPCAUtils method)

[TestCli](#) (class in moseq2\_pca.tests.integration\_tests.test\_cli)

[TestGUI](#) (class in moseq2\_pca.tests.integration\_tests.test\_gui)

[TestPCAUtils](#) (class in moseq2\_pca.tests.unit\_tests.test\_pca\_util)

[TestUtils](#) (class in moseq2\_pca.tests.unit\_tests.test\_util)

[TestViz](#) (class in moseq2\_pca.tests.unit\_tests.test\_viz)

[train\\_pca\\_command\(\)](#) (in module moseq2\_pca.gui)

[train\\_pca\\_dask\(\)](#) (in module moseq2\_pca.pca.util)

[train\\_pca\\_wrapper\(\)](#) (in module moseq2\_pca.helpers.wrappers)



# Python Module Index

## *m*

[moseq2\\_pca](#)

[moseq2\\_pca.gui](#)

[moseq2\\_pca.helpers.data](#)

[moseq2\\_pca.helpers.wrappers](#)

[moseq2\\_pca.pca.util](#)

[moseq2\\_pca.tests.integration\\_tests.test\\_cli](#)

[moseq2\\_pca.tests.integration\\_tests.test\\_gui](#)

[moseq2\\_pca.tests.unit\\_tests.test\\_pca\\_util](#)

[moseq2\\_pca.tests.unit\\_tests.test\\_util](#)

[moseq2\\_pca.tests.unit\\_tests.test\\_viz](#)

[moseq2\\_pca.util](#)

[moseq2\\_pca.viz](#)