

Python Documentation

version

April 14, 2020

Contents

Welcome to moseq2-pca's documentation!	1
moseq2-pca	1
moseq2_pca package	1
Subpackages	1
moseq2_pca.helpers package	1
Submodules	1
moseq2_pca.helpers.data module	1
moseq2_pca.helpers.wrappers module	1
Module contents	2
moseq2_pca.pca package	2
Submodules	2
moseq2_pca.pca.util module	2
Module contents	5
moseq2_pca.tests package	5
Subpackages	5
moseq2_pca.tests.integration_tests package	5
Submodules	5
moseq2_pca.tests.integration_tests.test_cli module	5
moseq2_pca.tests.integration_tests.test_gui module	5
Module contents	5
moseq2_pca.tests.unit_tests package	5
Submodules	5
moseq2_pca.tests.unit_tests.test_pca_util module	5
moseq2_pca.tests.unit_tests.test_util module	5
moseq2_pca.tests.unit_tests.test_viz module	6
Module contents	6
Module contents	6
Submodules	6
moseq2_pca.cli module	6
moseq2_pca.gui module	6
moseq2_pca.util module	7
moseq2_pca.viz module	10
Module contents	10
Indices and tables	10
Index	13
Python Module Index	17

Welcome to moseq2-pca's documentation!

moseq2-pca

moseq2_pca package

Subpackages

moseq2_pca.helpers package

Submodules

moseq2_pca.helpers.data module

`moseq2_pca.helpers.data.get_pca_yaml_data` (pca_yaml)

Reads PCA yaml file and returns metadata

Parameters: `pca_yaml` (**str**) (path to pca.yaml)

Returns: `use_fft` (**bool**) (indicates whether to use FFT) `clean_params` (**dict**) (dict of image filtering parameters) `mask_params` (**dict**) (dict of mask parameters) `missing_data` (**bool**) (indicates whether to use mask_params)

`moseq2_pca.helpers.data.load_pcs_for_cp` (pca_file_components, config_data)

Load computed Principal Components.

Parameters:

- `pca_file_components` (**str**) (path to pca h5 file to read PCs)
- `config_data` (**dict**) (config parameters)

Returns: `pca_components` (**str**) (path to pca components) `changepoint_params` (**dict**) (dict of relevant changepoint parameters) `cluster` (**dask Cluster**) (Dask Cluster object.) `client` (**dask Client**) (Dask Client Object) `missing_data` (**bool**) (Indicates whether to use mask_params) `mask_params` (**dict**) (Mask parameters to use when computing CPs)

`moseq2_pca.helpers.data.setup_cp_command` (input_dir, config_data, output_dir, output_file, output_directory)

Helper function for changepoints_wrapper to perform data-path existence checks.

Parameters:

- `input_dir` (**int**) (path to directory containing all h5+yaml files)
- `config_data` (**dict**) (dict of relevant PCA parameters (image filtering etc.))
- `output_dir` (**str**) (path to directory to store PCA data)
- `output_file` (**str**) (pca model filename)
- `output_directory` (**str**) (alternative output_dir)

Returns: `config_data` (**dict**) (updated config_data dict with the proper paths) `pca_file_components` (**str**) (path to trained pca file) `pca_file_scores` (**str**) (path to pca_scores file) `h5s` (**list**) (list of relevant pca h5 files) `yamls` (**list**) (list of relevant pca metadata yaml files) `save_file` (**str**) (path to save changepoints)

moseq2_pca.helpers.wrappers module

`moseq2_pca.helpers.wrappers.apply_pca_wrapper` (input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)

Wrapper function to obtain PCA Scores.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.compute_changepoints_wrapper` (input_dir, config_data, output_dir, output_file, gui=False, output_directory=None)

Wrapper function to compute model-free (PCA based) Changepoints.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.train_pca_wrapper` (input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)

Wrapper function to train PCA.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

Module contents

`moseq2_pca.pca` package

Submodules

`moseq2_pca.pca.util` module

`moseq2_pca.pca.util.apply_pca_dask` (pca_components, h5s, yamls, use_fft, clean_params, save_file, chunk_size, mask_params, missing_data, client, fps=30, gui=False)

“Apply” trained PCA on input frame data to obtain PCA Scores using Distributed Dask cluster.

Parameters:

- **pca_components (np.array)** (*array of computed Principal Components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **use_fft (bool)** (*indicate whether to use 2D-FFT*)
- **clean_params (dict)** (*dictionary containing filtering options*)
- **save_file (str)** (*path to pca_scores filename to save*)
- **chunk_size (int)** (*size of chunks to process*)
- **mask_params (dict)** (*dictionary of masking parameters (if missing data)*)
- **missing_data (bool)** (*indicates whether to use mask arrays.*)
- **fps (int)** (*frames per second*)

Returns:

Return type: None

`moseq2_pca.pca.util.apply_pca_local` (pca_components, h5s, yamls, use_fft, clean_params, save_file, chunk_size, mask_params, missing_data, fps=30)

“Apply” trained PCA on input frame data to obtain PCA Scores using local cluster/platform.

Parameters:

- **pca_components (np.array)** (*array of computed Principal Components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **use_fft (bool)** (*indicate whether to use 2D-FFT*)
- **clean_params (dict)** (*dictionary containing filtering options*)
- **save_file (str)** (*path to pca_scores filename to save*)
- **chunk_size (int)** (*size of chunks to process*)
- **mask_params (dict)** (*dictionary of masking parameters (if missing data)*)
- **missing_data (bool)** (*indicates whether to use mask arrays.*)
- **fps (int)** (*frames per second*)

Returns:

Return type: None

`moseq2_pca.pca.util.get_changepoints_dask` (changepoint_params, pca_components, h5s, yamls, save_file, chunk_size, mask_params, missing_data, client, fps=30, pca_scores=None, progress_bar=False, gui=False)

Computes model-free changepoints using PCs and PC Scores on distributed dask cluster.

Parameters:

- **changepoint_params (dict)** (*dict of changepoint parameters*)
- **pca_components (np.array)** (*computed principal components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **save_file (str)** (*path to save changepoint files*)
- **chunk_size (int)** (*size of chunks to process in dask.*)
- **mask_params (dict)** (*dict of missing_data mask parameters.*)
- **missing_data (bool)** (*indicate whether to use mask_params*)
- **client (dask Client)** (*initialized Dask Client object*)
- **fps (int)** (*frames per second*)
- **pca_scores (np.array)** (*computed principal component scores*)
- **progress_bar (bool)** (*display progress bar*)
- **gui (bool)** (*indicate GUI use*)

Returns:

Return type: None

`moseq2_pca.pca.util.mask_data` (`original_data`, `mask`, `new_data`)
Create a mask subregion given a boolean mask if missing data flag is used.

Parameters:

- **original_data (3d np.ndarray)** (*input frames*)
- **mask (3d boolean np.ndarray)** (*mask array*)
- **new_data (3d np.ndarray)** (*frames to use*)

Returns: output (3d np.ndarray)

Return type: masked data array

`moseq2_pca.pca.util.train_pca_dask` (`dask_array`, `clean_params`, `use_fft`, `rank`, `cluster_type`, `client`, `workers`, `cache`, `mask=None`, `iters=10`, `recon_pcs=10`, `min_height=10`, `max_height=100`)
Train PCA using dask arrays.

Parameters:

- **dask_array (dask array)** (*chunked frames to train PCA*)
- **clean_params (dict)** (*dictionary containing filtering parameters*)
- **use_fft (bool)** (*indicates whether to use 2d-FFT on images.*)
- **rank (int)** (*Matrix rank to use*)
- **cluster_type (str)** (*indicates which cluster to use.*)
- **client (Dask.Client)** (*client object to execute dask operations*)
- **workers (int)** (*number of dask workers*)
- **cache (str)** (*path to cache directory*)
- **mask (dask array)** (*dask array of masked data if missing_data parameter==True*)
- **iters (int)** (*number of SVD iterations*)
- **recon_pcs (int)** (*number of PCs to reconstruct. (if missing_data = True)*)
- **min_height (int)** (*minimum mouse height from floor in (mm)*)
- **max_height (int)** (*maximum mouse height from floor in (mm)*)

Returns: output_dict (dict)

Return type: dictionary containing PCA training results.

Module contents

moseq2_pca.tests package

Subpackages

moseq2_pca.tests.integration_tests package

Submodules

moseq2_pca.tests.integration_tests.test_cli module

```
class moseq2_pca.tests.integration_tests.test_cli.TestCli (methodName='runTest')
  Bases: unittest.case.TestCase

  test_apply_pca ()

  test_clip_scores ()

  test_compute_changepoints ()

  test_train_pca ()
```

moseq2_pca.tests.integration_tests.test_gui module

```
class moseq2_pca.tests.integration_tests.test_gui.TestGUI (methodName='runTest')
  Bases: unittest.case.TestCase

  test_apply_pca_command ()

  test_compute_changepoints_command ()

  test_train_pca_command ()
```

Module contents

moseq2_pca.tests.unit_tests package

Submodules

moseq2_pca.tests.unit_tests.test_pca_util module

```
class moseq2_pca.tests.unit_tests.test_pca_util.TestPCAUtils (methodName='runTest')
  Bases: unittest.case.TestCase

  test_mask_data ()

  test_train_pca_dask ()
```

moseq2_pca.tests.unit_tests.test_util module

```
class moseq2_pca.tests.unit_tests.test_util.TestUtils (methodName='runTest')
  Bases: unittest.case.TestCase

  test_clean_frames ()
```

Welcome to moseq2-pca's documentation!

```
test_gauss_smooth ()
test_gaussian_kernel1d ()
test_get_changepoints ()
test_get_metadata_path ()
test_get_rps ()
test_get_rsp_dask ()
test_get_timestamp_path ()
test_initialize_dask ()
test_insert_nans ()
test_read_yaml ()
test_recursive_find_h5s ()
test_select_strel ()
```

moseq2_pca.tests.unit_tests.test_viz module

```
class moseq2_pca.tests.unit_tests.test_viz.TestViz (methodName='runTest')
  Bases: unittest.case.TestCase

  changepoint_dist ()

  test_display_components ()
    cmap = 'gray' im_size = int(np.sqrt(components.shape[1])) plotv = components.reshape((-1, im_size, im_size))
    plotv = skimage.util.montage(plotv)
    plt.switch_backend('agg')
    fig, ax = plt.subplots(1, 1, figsize=(10, 10)) plt.imshow(plotv, cmap=cmap) plt.xticks([]) plt.yticks([])

  test_scee_plot ()
```

Module contents

Module contents

Submodules

moseq2_pca.cli module

```
moseq2_pca.cli.new_init (self, *args, **kwargs)
```

moseq2_pca.gui module

```
moseq2_pca.gui.apply_pca_command (input_dir, index_file, config_file, output_dir, output_file,
output_directory=None)
```

Compute PCA Scores given trained PCA using Jupyter Notebook.

Parameters:

- **input_dir (str)** (*path to directory containing training data*)
- **index_file (str)** (*path to index file.*)
- **config_file (str)** (*path to config file*)
- **output_dir (str)** (*path to output pca directory*)
- **output_file (str)** (*name of output pca file.*)
- **output_directory (str)** (*alternative output directory path*)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.compute_changepoints_command` (input_dir, config_file, output_dir, output_file, output_directory=None)

Compute Changepoint distribution using Jupyter Notebook.

Parameters:

- **input_dir (str)** (*path to directory containing training data*)
- **config_file (str)** (*path to config file*)
- **output_dir (str)** (*path to output pca directory*)
- **output_file (str)** (*name of output pca file.*)
- **output_directory (str)** (*alternative output directory path*)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.train_pca_command` (input_dir, config_file, output_dir, output_file, output_directory=None)

Train PCA through Jupyter notebook, and updates config file.

Parameters:

- **input_dir (str)** (*path to directory containing training data*)
- **config_file (str)** (*path to config file*)
- **output_dir (str)** (*path to output pca directory*)
- **output_file (str)** (*name of output pca file.*)
- **output_directory (str)** (*alternative output directory path*)

Returns:

Return type: None

moseq2_pca.util module

`moseq2_pca.util.clean_frames` (frames, medfilter_space=None, gaussfilter_space=None, medfilter_time=None, gaussfilter_time=None, detrend_time=None, tailfilter=None, tail_threshold=5)

Filters spatial/temporal noise from frames using Median and Gaussian filters, given kernel sizes for each respective requested filter.

Parameters:

- **frames (3D numpy array)** (*frames to filter.*)
- **medfilter_space (list)** (*median spatial filter kernel.*)
- **gaussfilter_space (list)** (*gaussian spatial filter kernel.*)
- **medfilter_time (list)** (*median temporal filter.*)
- **gaussfilter_time (list)** (*gaussian temporal filter.*)
- **detrend_time (int)** (*number of frames to lag for.*)
- **tailfilter (int)** (*size of tail-filter kernel.*)
- **tail_threshold (int)** (*threshold value to use for tail filtering*)

Returns: out (3D numpy array)

Return type: filtered frames.

`moseq2_pca.util.command_with_config` (config_file_param_name)

`moseq2_pca.util.gauss_smooth` (signal, win_length=None, sig=1.5, kernel=None)
Perform Gaussian Smoothing on a 1D signal.

Parameters:

- **signal (1d numpy array)** (*signal to perform smoothing*)
- **win_length (int)** (*window_size for gaussian kernel filter*)
- **sig (float)** (*variance of 1d gaussian kernel.*)
- **kernel (tuple)** (*kernel size to use for smoothing*)

Returns: **result (1d numpy array)**

Return type: smoothed signal

`moseq2_pca.util.gaussian_kernel1d` (n=None, sig=3)
Get 1D gaussian kernel.

Parameters:

- **n (int)** (*number of points to use.*)
- **sig (int)** (*variance of kernel to use.*)

Returns: **kernel (1d array)**

Return type: 1D numpy kernel.

`moseq2_pca.util.get_changepoints` (scores, k=5, sigma=3, peak_height=0.5, peak_neighbors=1, baseline=True, timestamps=None)
Compute changepoints distribution and CP Curve.

Parameters:

- **scores (3D numpy array)** (*nframes * r * c*)
- **k (int)** (*klags - Lag to use for derivative calculation.*)
- **sigma (int)** (*Standard deviation of gaussian smoothing filter.*)
- **peak_height (float)** (*user-defined peak Changepoint length.*)
- **peak_neighbors (int)** (*number of peaks in the CP curve.*)
- **baseline (bool)** (*normalize data.*)
- **timestamps (array)** (*loaded timestamps.*)

Returns: **cps (numpy array)** (*array of values for CP curve*) **normed_df (numpy array)** (*array of values for bar plot*)

`moseq2_pca.util.get_metadata_path` (h5file)
Return path within h5 file that contains the kinect extraction metadata.

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to acquisition metadata within h5 file.

`moseq2_pca.util.get_rps` (frames, rps=600, normalize=True)
Get random projections of frames.

Parameters:

- **frames (2D or 3D numpy array)** (*Frames to get dimensions from.*)
- **rps (int)** (*Number of random projections.*)
- **normalize (bool)** (*indicates whether to normalize frames.*)

Returns: **rproj (2D or 3D numpy array)**

Return type: Computed random projections with same shape as frames

`moseq2_pca.util.get_timestamp_path` (h5file)
Return path within h5 file that contains the kinect timestamps

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to metadata timestamps within h5 file

`moseq2_pca.util.initialize_dask` (nworkers=50, processes=1, memory='4GB', cores=1, wall_time='01:00:00', queue='debug', local_processes=False, cluster_type='local', scheduler='distributed', timeout=10, cache_path='/Users/aymanzeine/moseq2_pca', **kwargs)
Initialize dask client, cluster, workers, etc.

Parameters:

- **nworkers (int)** (*number of dask workers to initialize*)
- **processes (int)** (*number of processes per worker*)
- **memory (str)** (*amount of memory to allocate to dask cluster*)
- **cores (int)** (*number of cores to use.*)
- **wall_time (str)** (*amount of time to allow program to run*)
- **queue (str)** (*logging mode*)
- **local_processes (bool)** (*indicate whether the processes are local*)
- **cluster_type (str)** (*indicate what cluster to use*)
- **scheduler (str)** (*indicate what scheduler to use*)
- **timeout (int)** (*number of worker timeouts to allow*)
- **cache_path (str or Pathlike)** (*path to store cached data*)
- **kwargs** (*extra keyword arguments*)

Returns: **client (dask Client)** (*initialized Client*) **cluster (dask Cluster)** (*initialized Cluster*) **workers (dask Workers)** (*intialized workers*) **cache (dask Chest)** (*initialized Chest (cache) object*)

`moseq2_pca.util.insert_nans` (timestamps, data, fps=30)
Fills NaN values with 0 in timestamps.

Parameters:

- **timestamps (1D array)** (*timestamp time-strs*)
- **data (1D array)** (*timestamp values*)
- **fps (int)** (*frames per second*)

Returns: **filled_data (1D array)** (*filled missing timestamp values.*) **data_idx (1D array)** (*indices of inserted 0s*) **filled_timestamps (1D array)** (*filled timestamp-strs*)

`moseq2_pca.util.read_yaml` (yaml_file)
Reads yaml file and returns dictionary representation of file contents.

Parameters: **yaml_file (str)** (*path to yaml file*)

Returns: **return_dict (dict)**

Return type: dict of yaml file contents

`moseq2_pca.util.recursive_find_h5s` (root_dir='/Users/aymanzeine/Desktop/moseq/moseq2-pca/docs', ext='.h5', yaml_string='{}.yaml')
Recursively find h5 files, along with yaml files with the same basename

Parameters:

- **root_dir (str or os.Pathlike)** (*path to directory to start recursive search*)
- **ext (str)** (*extension to search for, e.g. .h5*)
- **yaml_string (str)** (*a format to use to name yaml files*)

Returns: **h5s (list)** (*list of h5 file paths*) **dicts (list)** (*list of metadata file paths*) **yamls (list)** (*list of yaml file paths*)

`moseq2_pca.util.recursively_load_dict_contents_from_group` (h5file, path)
Reads all contents from h5 and returns them in a nested dict object.

Parameters:

- **h5file (str)** (*path to h5 file*)
- **path (str)** (*path to group within h5 file*)

Returns: **ans (dict)****Return type:** dictionary of all h5 group contents

`moseq2_pca.util.select_strel` (string='e', size=10, 10)
 Selects Structuring Element Shape

Parameters:

- **string (str)** (*e for Ellipse, r for Rectangle*)
- **size (tuple)** (*size of StructuringElement*)

Returns: **strel (cv2.StructuringElement)****Return type:** returned StructuringElement with specified size.

`moseq2_pca.util.shutdown_dask` (scheduler)

Graceful shutdown dask scheduler. source:
<https://github.com/dask/distributed/issues/1703#issuecomment-361291492>

Parameters: **scheduler (dask Scheduler)** (*scheduler to shutdown.*)**Returns:****Return type:** None***moseq2_pca.viz module***

`moseq2_pca.viz.changepoint_dist` (cps, headless=False)
 Creates bar plot describing computed Changepoint Distribution.

Parameters:

- **cps (np.ndarray)** (*changepoints to graph*)
- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.display_components` (components, cmap='gray', headless=False)
 Creates grid of computed Principal Components.

Parameters:

- **components (np.ndarray)** (*components to graph*)
- **cmap (str)** (*color map to use*)
- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.scree_plot` (explained_variance_ratio, headless=False)
 Creates Scree plot describing principal components.

Parameters:

- **explained_variance_ratio (np.array)** (*explained variance ratio of each principal component*)
- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)****Return type:** figure to save/graph***Module contents*****Indices and tables**

- **genindex**
- **modindex**

- **search**

Index

A

[apply_pca_command\(\)](#) (in module [moseq2_pca.gui](#))
[apply_pca_dask\(\)](#) (in module [moseq2_pca.pca.util](#))
[apply_pca_local\(\)](#) (in module [moseq2_pca.pca.util](#))
[apply_pca_wrapper\(\)](#) (in module [moseq2_pca.helpers.wrappers](#))

C

[changepoint_dist\(\)](#) (in module [moseq2_pca.viz](#))
 ([moseq2_pca.tests.unit_tests.test_viz.TestViz](#) method)
[clean_frames\(\)](#) (in module [moseq2_pca.util](#))
[command_with_config\(\)](#) (in module [moseq2_pca.util](#))
[compute_changepoints_command\(\)](#) (in module [moseq2_pca.gui](#))
[compute_changepoints_wrapper\(\)](#) (in module [moseq2_pca.helpers.wrappers](#))

D

[display_components\(\)](#) (in module [moseq2_pca.viz](#))

G

[gauss_smooth\(\)](#) (in module [moseq2_pca.util](#))
[gaussian_kernel1d\(\)](#) (in module [moseq2_pca.util](#))
[get_changepoints\(\)](#) (in module [moseq2_pca.util](#))
[get_changepoints_dask\(\)](#) (in module [moseq2_pca.pca.util](#))
[get_metadata_path\(\)](#) (in module [moseq2_pca.util](#))
[get_pca_yaml_data\(\)](#) (in module [moseq2_pca.helpers.data](#))
[get_rps\(\)](#) (in module [moseq2_pca.util](#))
[get_timestamp_path\(\)](#) (in module [moseq2_pca.util](#))

I

[initialize_dask\(\)](#) (in module [moseq2_pca.util](#))
[insert_nans\(\)](#) (in module [moseq2_pca.util](#))

L

[load_pcs_for_cp\(\)](#) (in module [moseq2_pca.helpers.data](#))

M

[mask_data\(\)](#) (in module [moseq2_pca.pca.util](#))
module

[moseq2_pca](#)
[moseq2_pca.cli](#)
[moseq2_pca.gui](#)
[moseq2_pca.helpers](#)
[moseq2_pca.helpers.data](#)
[moseq2_pca.helpers.wrappers](#)
[moseq2_pca.pca](#)
[moseq2_pca.pca.util](#)
[moseq2_pca.tests](#)
[moseq2_pca.tests.integration_tests](#)
[moseq2_pca.tests.integration_tests.test_cli](#)
[moseq2_pca.tests.integration_tests.test_gui](#)
[moseq2_pca.tests.unit_tests](#)
[moseq2_pca.tests.unit_tests.test_pca_util](#)
[moseq2_pca.tests.unit_tests.test_util](#)
[moseq2_pca.tests.unit_tests.test_viz](#)
[moseq2_pca.util](#)
[moseq2_pca.viz](#)

moseq2_pca

[module](#)

moseq2_pca.cli

[module](#)

moseq2_pca.gui

[module](#)

moseq2_pca.helpers

[module](#)

moseq2_pca.helpers.data

[module](#)

moseq2_pca.helpers.wrappers

[module](#)

moseq2_pca.pca

[module](#)

moseq2_pca.pca.util

[module](#)

moseq2_pca.tests

[module](#)

moseq2_pca.tests.integration_tests

[module](#)

moseq2_pca.tests.integration_tests.test_cli

[module](#)

moseq2_pca.tests.integration_tests.test_gui

[module](#)

moseq2_pca.tests.unit_tests

[module](#)

moseq2_pca.tests.unit_tests.test_pca_util

[module](#)

moseq2_pca.tests.unit_tests.test_util

module

moseq2_pca.tests.unit_tests.test_viz

module

moseq2_pca.util

module

moseq2_pca.viz

module

N

`new_init()` (in module `moseq2_pca.cli`)

R

`read_yaml()` (in module `moseq2_pca.util`)

`recursive_find_h5s()` (in module `moseq2_pca.util`)

`recursively_load_dict_contents_from_group()` (in module `moseq2_pca.util`)

S

`scree_plot()` (in module `moseq2_pca.viz`)

`select_strel()` (in module `moseq2_pca.util`)

`setup_cp_command()` (in module `moseq2_pca.helpers.data`)

`shutdown_dask()` (in module `moseq2_pca.util`)

T

`test_apply_pca()`
(`moseq2_pca.tests.integration_tests.test_cli.TestCli` method)

`test_apply_pca_command()`
(`moseq2_pca.tests.integration_tests.test_gui.TestGUI` method)

`test_clean_frames()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_clip_scores()`
(`moseq2_pca.tests.integration_tests.test_cli.TestCli` method)

`test_compute_changepoints()`
(`moseq2_pca.tests.integration_tests.test_cli.TestCli` method)

`test_compute_changepoints_command()`
(`moseq2_pca.tests.integration_tests.test_gui.TestGUI` method)

`test_display_components()`
(`moseq2_pca.tests.unit_tests.test_viz.TestViz` method)

`test_gauss_smooth()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_gaussian_kernel1d()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_get_changepoints()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_get_metadata_path()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_get_rps()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_get_rsp_dask()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_get_timestamp_path()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_initialize_dask()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_insert_nans()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_mask_data()` (`moseq2_pca.tests.unit_tests.test_pca_util.TestPCAUtils` method)

`test_read_yaml()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_recursive_find_h5s()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_scree_plot()`
(`moseq2_pca.tests.unit_tests.test_viz.TestViz` method)

`test_select_strel()`
(`moseq2_pca.tests.unit_tests.test_util.TestUtils` method)

`test_train_pca()`
(`moseq2_pca.tests.integration_tests.test_cli.TestCli` method)

`test_train_pca_command()`
(`moseq2_pca.tests.integration_tests.test_gui.TestGUI` method)

`test_train_pca_dask()` (`moseq2_pca.tests.unit_tests.test_pca_util.TestPCAUtils` method)

`TestCli` (class in `moseq2_pca.tests.integration_tests.test_cli`)

`TestGUI` (class in `moseq2_pca.tests.integration_tests.test_gui`)

`TestPCAUtils` (class in `moseq2_pca.tests.unit_tests.test_pca_util`)

`TestUtils` (class in `moseq2_pca.tests.unit_tests.test_util`)

TestViz (class in moseq2_pca.tests.unit_tests.test_viz)
train_pca_command() (in module moseq2_pca.gui)
train_pca_dask() (in module moseq2_pca.pca.util)
train_pca_wrapper() (in module moseq2_pca.helpers.wrappers)

Python Module Index

m

[moseq2_pca](#)

[moseq2_pca.cli](#)

[moseq2_pca.gui](#)

[moseq2_pca.helpers](#)

[moseq2_pca.helpers.data](#)

[moseq2_pca.helpers.wrappers](#)

[moseq2_pca.pca](#)

[moseq2_pca.pca.util](#)

[moseq2_pca.tests](#)

[moseq2_pca.tests.integration_tests](#)

[moseq2_pca.tests.integration_tests.test_cli](#)

[moseq2_pca.tests.integration_tests.test_gui](#)

[moseq2_pca.tests.unit_tests](#)

[moseq2_pca.tests.unit_tests.test_pca_util](#)

[moseq2_pca.tests.unit_tests.test_util](#)

[moseq2_pca.tests.unit_tests.test_viz](#)

[moseq2_pca.util](#)

[moseq2_pca.viz](#)