

Python Documentation

version

May 13, 2020

Contents

Welcome to moseq2-pca's documentation!	1
moseq2_pca package	1
CLI Module	1
moseq2-pca	1
apply-pca	1
clip-scores	2
compute-changepoints	2
train-pca	3
GUI Module	5
Utilities Module	5
Visualization Module	8
Subpackages	9
moseq2_pca.helpers package	9
Helpers - Data Module	9
Helpers - Wrapper Module	9
moseq2_pca.pca package	10
PCA - Utilities Module	10
Indices and tables	13
Index	15
Python Module Index	21

Welcome to moseq2-pca's documentation!

moseq2_pca package

CLI Module

moseq2-pca

```
moseq2-pca [OPTIONS] COMMAND [ARGS]...
```

apply-pca

```
moseq2-pca apply-pca [OPTIONS]
```

Options

- i, --input-dir** <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]
- cluster-type** <cluster_type>
Cluster type [default: local]
Options: local|slurm|nodask
- o, --output-dir** <output_dir>
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca]
- output-file** <output_file>
Name of h5 file for storing pca results [default: pca_scores]
- h5-path** <h5_path>
Path to data in h5 files [default: /frames]
- h5-mask-path** <h5_mask_path>
Path to log-likelihood mask in h5 files [default: /frames_mask]
- pca-path** <pca_path>
Path to pca components [default: /components]
- pca-file** <pca_file>
Path to PCA results
- chunk-size** <chunk_size>
Number of frames per chunk [default: 4000]
- fill-gaps** <fill_gaps>
Fill dropped frames with nans [default: True]
- fps** <fps>
Fps (only used if no timestamps found) [default: 30]
- detrend-window** <detrend_window>
Length of detrend window (in seconds, 0 for no detrending) [default: 0]
- config-file** <config_file>
Path to configuration file
- d, --dask-cache-path** <dask_cache_path>
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]
- q, --queue** <queue>
Cluster queue/partition for submitting jobs [default: debug]
- n, --nworkers** <nworkers>
Number of workers [default: 10]

Welcome to moseq2-pca's documentation!

-c, --cores <cores>
Number of cores per worker [default: 1]

-p, --processes <processes>
Number of processes to run on each worker [default: 1]

-m, --memory <memory>
RAM usage per workers [default: 15GB]

-w, --wall-time <wall_time>
Wall time for workers [default: 06:00:00]

--timeout <timeout>
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

clip-scores

Clips PCA scores from the beginning or end

Args:

pca_file (string): Path to PCA scores **clip_samples** (int): number of samples to clip from beginning or end
from_end (bool): if true clip from end rather than beginning

Note that scores are modified *in place*.

```
moseq2-pca clip-scores [OPTIONS] PCA_FILE CLIP_SAMPLES
```

Options

--from-end
[default: False]

Arguments

PCA_FILE
Required argument

CLIP_SAMPLES
Required argument

compute-changepoints

```
moseq2-pca compute-changepoints [OPTIONS]
```

Options

-i, --input-dir <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

-o, --output-dir <output_dir>
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca/]

--output-file <output_file>
Name of h5 file for storing pca results [default: changepoints]

--cluster-type <cluster_type>
Cluster type [default: local]
Options: local|slurm

--pca-file-components <pca_file_components>
Path to PCA components

--pca-file-scores <pca_file_scores>
Path to PCA results

--pca-path <pca_path>
Path to pca components [default: /components]

--neighbors <neighbors>

Welcome to moseq2-pca's documentation!

Neighbors to use for peak identification [default: 1]

--threshold <threshold>
Peak threshold to use for changepoints [default: 0.5]

-k, --klags <klags>
Lag to use for derivative calculation [default: 6]

-s, --sigma <sigma>
Standard deviation of gaussian smoothing filter [default: 3.5]

-d, --dims <dims>
Number of random projections to use [default: 300]

--fps <fps>
Fps (only used if no timestamps found) [default: 30]

--h5-path <h5_path>
Path to data in h5 files [default: /frames]

--h5-mask-path <h5_mask_path>
Path to log-likelihood mask in h5 files [default: /frames_mask]

--chunk-size <chunk_size>
Number of frames per chunk [default: 4000]

--config-file <config_file>
Path to configuration file

--dask-cache-path <dask_cache_path>
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

--visualize-results <visualize_results>
Visualize results [default: True]

-q, --queue <queue>
Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
Number of workers [default: 10]

-c, --cores <cores>
Number of cores per worker [default: 1]

-p, --processes <processes>
Number of processes to run on each worker [default: 1]

-m, --memory <memory>
RAM usage per workers [default: 15GB]

-w, --wall-time <wall_time>
Wall time for workers [default: 06:00:00]

--timeout <timeout>
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

train-pca

```
moseq2-pca train-pca [OPTIONS]
```

Options

-i, --input-dir <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

--cluster-type <cluster_type>
Cluster type [default: local]

Options: local|slurm

-o, --output-dir <output_dir>

Welcome to moseq2-pca's documentation!

Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca]

--gaussfilter-space <gaussfilter_space>
Spatial filter for data (Gaussian) [default: 1.5, 1]

--gaussfilter-time <gaussfilter_time>
Temporal filter for data (Gaussian) [default: 0]

--medfilter-space <medfilter_space>
Median spatial filter [default: 0]

--medfilter-time <medfilter_time>
Median temporal filter [default: 0]

--missing-data
Use missing data PCA [default: False]

--missing-data-iters <missing_data_iters>
Missing data PCA iterations [default: 10]

--mask-threshold <mask_threshold>
Threshold for mask (missing data only) [default: -16]

--mask-height-threshold <mask_height_threshold>
Threshold for mask based on floor height [default: 5]

--min-height <min_height>
Min mouse height from floor (mm) [default: 10]

--max-height <max_height>
Max mouse height from floor (mm) [default: 100]

--tailfilter-size <tailfilter_size>
Tail filter size [default: 9, 9]

--tailfilter-shape <tailfilter_shape>
Tail filter shape [default: ellipse]

--use-fft
Use 2D fft [default: False]

--recon-pcs <recon_pcs>
Number of PCs to use for missing data reconstruction [default: 10]

--rank <rank>
Rank for compressed SVD (generally>>nPCS) [default: 50]

--output-file <output_file>
Name of h5 file for storing pca results [default: pca]

--chunk-size <chunk_size>
Number of frames per chunk [default: 4000]

--visualize-results <visualize_results>
Visualize results [default: True]

--config-file <config_file>
Path to configuration file

-d, --dask-cache-path <dask_cache_path>
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

--local-processes <local_processes>
Use processes with local scheduler [default: True]

-q, --queue <queue>
Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
Number of workers [default: 10]

-c, --cores <cores>
Number of cores per worker [default: 1]

Welcome to moseq2-pca's documentation!

-p, --processes <processes>
Number of processes to run on each worker [default: 1]

-m, --memory <memory>
Total RAM usage per worker [default: 15GB]

-w, --wall-time <wall_time>
Wall time for workers [default: 06:00:00]

--timeout <timeout>
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

GUI Module

`moseq2_pca.gui.apply_pca_command(input_dir, index_file, config_file, output_dir, output_file, output_directory=None)`

Compute PCA Scores given trained PCA using Jupyter Notebook.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **index_file (str)** (path to index file.)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.compute_changepoints_command(input_dir, config_file, output_dir, output_file, output_directory=None)`

Compute Changepoint distribution using Jupyter Notebook.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.train_pca_command(input_dir, config_file, output_dir, output_file, output_directory=None)`

Train PCA through Jupyter notebook, and updates config file.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns:

Return type: None

Utilities Module

`moseq2_pca.util.clean_frames` (frames, medfilter_space=None, gaussfilter_space=None, medfilter_time=None, gaussfilter_time=None, detrend_time=None, tailfilter=None, tail_threshold=5)

Filters spatial/temporal noise from frames using Median and Gaussian filters, given kernel sizes for each respective requested filter.

Parameters:

- **frames (3D numpy array)** (*frames to filter.*)
- **medfilter_space (list)** (*median spatial filter kernel.*)
- **gaussfilter_space (list)** (*gaussian spatial filter kernel.*)
- **medfilter_time (list)** (*median temporal filter.*)
- **gaussfilter_time (list)** (*gaussian temporal filter.*)
- **detrend_time (int)** (*number of frames to lag for.*)
- **tailfilter (int)** (*size of tail-filter kernel.*)
- **tail_threshold (int)** (*threshold value to use for tail filtering*)

Returns: out (3D numpy array)

Return type: filtered frames.

`moseq2_pca.util.command_with_config` (config_file_param_name)

`moseq2_pca.util.gauss_smooth` (signal, win_length=None, sig=1.5, kernel=None)
Perform Gaussian Smoothing on a 1D signal.

Parameters:

- **signal (1d numpy array)** (*signal to perform smoothing*)
- **win_length (int)** (*window_size for gaussian kernel filter*)
- **sig (float)** (*variance of 1d gaussian kernel.*)
- **kernel (tuple)** (*kernel size to use for smoothing*)

Returns: result (1d numpy array)

Return type: smoothed signal

`moseq2_pca.util.gaussian_kernel1d` (n=None, sig=3)
Get 1D gaussian kernel.

Parameters:

- **n (int)** (*number of points to use.*)
- **sig (int)** (*variance of kernel to use.*)

Returns: kernel (1d array)

Return type: 1D numpy kernel.

`moseq2_pca.util.get_changepoints` (scores, k=5, sigma=3, peak_height=0.5, peak_neighbors=1, baseline=True, timestamps=None)

Compute changepoints distribution and CP Curve.

Parameters:

- **scores (3D numpy array)** (*nframes * r * c*)
- **k (int)** (*klags - Lag to use for derivative calculation.*)
- **sigma (int)** (*Standard deviation of gaussian smoothing filter.*)
- **peak_height (float)** (*user-defined peak Changepoint length.*)
- **peak_neighbors (int)** (*number of peaks in the CP curve.*)
- **baseline (bool)** (*normalize data.*)
- **timestamps (array)** (*loaded timestamps.*)

Returns: cps (2D numpy array) (array of values for CP curve) normed_df (1D numpy array) (array of values for bar plot)

`moseq2_pca.util.get_metadata_path` (h5file)

Return path within h5 file that contains the kinect extraction metadata.

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to acquisition metadata within h5 file.

`moseq2_pca.util.get_rps (frames, rps=600, normalize=True)`

Get random projections of frames.

Parameters:

- **frames (2D or 3D numpy array)** (*Frames to get dimensions from.*)
- **rps (int)** (*Number of random projections.*)
- **normalize (bool)** (*indicates whether to normalize frames.*)

Returns: **rproj (2D or 3D numpy array)**

Return type: Computed random projections with same shape as frames

`moseq2_pca.util.get_timestamp_path (h5file)`

Return path within h5 file that contains the kinect timestamps

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to metadata timestamps within h5 file

`moseq2_pca.util.initialize_dask (nworkers=50, processes=1, memory='4GB', cores=1, wall_time='01:00:00', queue='debug', local_processes=False, cluster_type='local', scheduler='distributed', timeout=10, cache_path='/Users/aymanzeine/moseq2_pca', **kwargs)`
Initialize dask client, cluster, workers, etc.

Parameters:

- **nworkers (int)** (*number of dask workers to initialize*)
- **processes (int)** (*number of processes per worker*)
- **memory (str)** (*amount of memory to allocate to dask cluster*)
- **cores (int)** (*number of cores to use.*)
- **wall_time (str)** (*amount of time to allow program to run*)
- **queue (str)** (*logging mode*)
- **local_processes (bool)** (*indicate whether the processes are local*)
- **cluster_type (str)** (*indicate what cluster to use*)
- **scheduler (str)** (*indicate what scheduler to use*)
- **timeout (int)** (*number of worker timeouts to allow*)
- **cache_path (str or Pathlike)** (*path to store cached data*)
- **kwargs** (*extra keyword arguments*)

Returns: **client (dask Client)** (*initialized Client*) **cluster (dask Cluster)** (*initialized Cluster*) **workers (dask Workers)** (*initialized workers or None if cluster_type = 'local'*) **cache (dask Chest)** (*initialized Chest (cache) object pointing to given cache path*)

`moseq2_pca.util.insert_nans (timestamps, data, fps=30)`

Fills NaN values with 0 in timestamps.

Parameters:

- **timestamps (1D array)** (*timestamp time-strs*)
- **data (1D array)** (*timestamp values*)
- **fps (int)** (*frames per second*)

Returns: **filled_data (1D array)** (*filled missing timestamp values.*) **data_idx (1D array)** (*indices of inserted 0s*) **filled_timestamps (1D array)** (*filled timestamp-strs*)

`moseq2_pca.util.read_yaml (yaml_file)`

Reads yaml file and returns dictionary representation of file contents.

Parameters: **yaml_file (str)** (*path to yaml file*)

Returns: **return_dict (dict)**

Return type: dict of yaml file contents

`moseq2_pca.util.recursive_find_h5s`

(`root_dir='/Users/aymanzeine/Desktop/moseq/moseq2-pca/docs'`, `ext='.h5'`,
`yaml_string='{ }.yaml'`)

Recursively find h5 files, along with yaml files with the same basename

Parameters:

- **root_dir (str or os.Pathlike)** (*path to directory to start recursive search*)

- **ext (str)** (*extension to search for, e.g. .h5*)

- **yaml_string (str)** (*a format to use to name yaml files*)

Returns: **h5s (list)** (*list of h5 file paths*) **dicts (list)** (*list of dicts containing metadata file contents*)

yamls (list) (*list of yaml file paths*)

`moseq2_pca.util.recursively_load_dict_contents_from_group` (`h5file`, `path`)

Reads all contents from h5 and returns them in a nested dict object.

Parameters:

- **h5file (str)** (*path to h5 file*)

- **path (str)** (*path to group within h5 file*)

Returns: **ans (dict)**

Return type: dictionary of all h5 group contents

`moseq2_pca.util.select_strel` (`string='e'`, `size=(10,10)`)

Selects Structuring Element Shape

Parameters:

- **string (str)** (*e for Ellipse, r for Rectangle*)

- **size (tuple)** (*size of StructuringElement*)

Returns: **strel (cv2.StructuringElement)**

Return type: returned StructuringElement with specified size.

`moseq2_pca.util.shutdown_dask` (`scheduler`)

Graceful shutdown dask scheduler. source:
<https://github.com/dask/distributed/issues/1703#issuecomment-361291492>

Parameters: **scheduler (dask Scheduler)** (*scheduler to shutdown.*)

Returns:

Return type: None

Visualization Module

`moseq2_pca.viz.changepoint_dist` (`cps`, `headless=False`)

Creates bar plot describing computed Changepoint Distribution.

Parameters:

- **cps (np.ndarray)** (*changepoints to graph*)

- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.display_components` (`components`, `cmap='gray'`, `headless=False`)

Creates grid of computed Principal Components.

Parameters:

- **components (2D np.ndarray)** (*components to graph*)

- **cmap (str)** (*color map to use*)

- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.scrree_plot` (*explained_variance_ratio, headless=False*)
Creates Scree plot describing principal components.

Parameters:

- **explained_variance_ratio (1D np.array)** (*explained variance ratio of each principal component*)
- **headless (bool)** (*trim first element in PC list*)

Returns: **plt (plt.figure)**

Return type: figure to save/graph

Subpackages

moseq2_pca.helpers package

Helpers - Data Module

`moseq2_pca.helpers.data.get_pca_yaml_data` (*pca_yaml*)
Reads PCA yaml file and returns metadata

Parameters: **pca_yaml (str)** (*path to pca.yaml*)

Returns: **use_fft (bool)** (*indicates whether to use FFT*) **clean_params (dict)** (*dict of image filtering parameters*) **mask_params (dict)** (*dict of mask parameters*) **missing_data (bool)** (*indicates whether to use mask_params*)

`moseq2_pca.helpers.data.load_pcs_for_cp` (*pca_file_components, config_data*)
Load computed Principal Components for Model-free Changepoint Analysis.

Parameters:

- **pca_file_components (str)** (*path to pca h5 file to read PCs*)
- **config_data (dict)** (*config parameters*)

Returns: **pca_components (str)** (*path to pca components*) **changepoint_params (dict)** (*dict of relevant changepoint parameters*) **cluster (dask Cluster)** (*Dask Cluster object.*) **client (dask Client)** (*Dask Client Object*) **missing_data (bool)** (*Indicates whether to use mask_params*) **mask_params (dict)** (*Mask parameters to use when computing CPs*)

`moseq2_pca.helpers.data.setup_cp_command` (*input_dir, config_data, output_dir, output_file, output_directory=None*)
Helper function for changepoints_wrapper to perform data-path existence checks.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)

Returns: **config_data (dict)** (*updated config_data dict with the proper paths*) **pca_file_components (str)** (*path to trained pca file*) **pca_file_scores (str)** (*path to pca_scores file*) **h5s (list)** (*list of relevant pca h5 files*) **yamls (list)** (*list of relevant pca metadata yaml files*) **save_file (str)** (*path to save changepoints*)

Helpers - Wrapper Module

`moseq2_pca.helpers.wrappers.apply_pca_wrapper` (*input_dir, config_data, output_dir, output_file, output_directory=None, gui=False*)
Wrapper function to obtain PCA Scores.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.compute_changepoints_wrapper` (input_dir, config_data, output_dir, output_file, gui=False, output_directory=None)

Wrapper function to compute model-free (PCA based) Changepoints.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.train_pca_wrapper` (input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)

Wrapper function to train PCA.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

moseq2_pca.pca package

PCA - Utilities Module

`moseq2_pca.pca.util.apply_pca_dask` (pca_components, h5s, yamls, use_fft, clean_params, save_file, chunk_size, mask_params, missing_data, client, fps=30, gui=False)

"Apply" trained PCA on input frame data to obtain PCA Scores using Distributed Dask cluster.

Parameters:

- **pca_components (np.array)** (*array of computed Principal Components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **use_fft (bool)** (*indicate whether to use 2D-FFT*)
- **clean_params (dict)** (*dictionary containing filtering options*)
- **save_file (str)** (*path to pca_scores filename to save*)
- **chunk_size (int)** (*size of chunks to process*)
- **mask_params (dict)** (*dictionary of masking parameters (if missing data)*)
- **missing_data (bool)** (*indicates whether to use mask arrays.*)
- **fps (int)** (*frames per second*)

Returns:

Return type: None

`moseq2_pca.pca.util.apply_pca_local` (`pca_components`, `h5s`, `yamls`, `use_fft`, `clean_params`, `save_file`, `chunk_size`, `mask_params`, `missing_data`, `fps=30`)

“Apply” trained PCA on input frame data to obtain PCA Scores using local cluster/platform.

Parameters:

- **pca_components (np.array)** (*array of computed Principal Components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **use_fft (bool)** (*indicate whether to use 2D-FFT*)
- **clean_params (dict)** (*dictionary containing filtering options*)
- **save_file (str)** (*path to pca_scores filename to save*)
- **chunk_size (int)** (*size of chunks to process*)
- **mask_params (dict)** (*dictionary of masking parameters (if missing data)*)
- **missing_data (bool)** (*indicates whether to use mask arrays.*)
- **fps (int)** (*frames per second*)

Returns:

Return type: None

`moseq2_pca.pca.util.get_changepoints_dask` (`changepoint_params`, `pca_components`, `h5s`, `yamls`, `save_file`, `chunk_size`, `mask_params`, `missing_data`, `client`, `fps=30`, `pca_scores=None`, `progress_bar=False`, `gui=False`)

Computes model-free changepoints using PCs and PC Scores on distributed dask cluster.

Parameters:

- **changepoint_params (dict)** (*dict of changepoint parameters*)
- **pca_components (np.array)** (*computed principal components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **save_file (str)** (*path to save changepoint files*)
- **chunk_size (int)** (*size of chunks to process in dask.*)
- **mask_params (dict)** (*dict of missing_data mask parameters.*)
- **missing_data (bool)** (*indicate whether to use mask_params*)
- **client (dask Client)** (*initialized Dask Client object*)
- **fps (int)** (*frames per second*)
- **pca_scores (np.array)** (*computed principal component scores*)
- **progress_bar (bool)** (*display progress bar*)
- **gui (bool)** (*indicate GUI use*)

Returns:

Return type: None

`moseq2_pca.pca.util.mask_data (original_data, mask, new_data)`
Create a mask subregion given a boolean mask if missing data flag is used.

Parameters:

- **original_data (3d np.ndarray)** (*input frames*)
- **mask (3d boolean np.ndarray)** (*mask array*)
- **new_data (3d np.ndarray)** (*frames to use*)

Returns: output (3d np.ndarray)

Return type: masked data array

`moseq2_pca.pca.util.train_pca_dask (dask_array, clean_params, use_fft, rank, cluster_type, client, workers, cache, mask=None, iters=10, recon_pcs=10, min_height=10, max_height=100)`
Train PCA using dask arrays.

Parameters:

- **dask_array (dask array)** (*chunked frames to train PCA*)
- **clean_params (dict)** (*dictionary containing filtering parameters*)
- **use_fft (bool)** (*indicates whether to use 2d-FFT on images.*)
- **rank (int)** (*Matrix rank to use*)
- **cluster_type (str)** (*indicates which cluster to use.*)
- **client (Dask.Client)** (*client object to execute dask operations*)
- **workers (int)** (*number of dask workers*)
- **cache (str)** (*path to cache directory*)
- **mask (dask array)** (*dask array of masked data if missing_data parameter==True*)
- **iters (int)** (*number of SVD iterations*)
- **recon_pcs (int)** (*number of PCs to reconstruct. (if missing_data = True)*)
- **min_height (int)** (*minimum mouse height from floor in (mm)*)
- **max_height (int)** (*maximum mouse height from floor in (mm)*)

Returns: output_dict (dict)

Return type: dictionary containing PCA training results.

Indices and tables

- `genindex`
- `modindex`
- `search`

Index

Symbols

		--h5-mask-path <h5_mask_path>	moseq2-pca-apply-pca	command line option
			moseq2-pca-compute-changepoints	command line option
--chunk-size <chunk_size>	moseq2-pca-apply-pca	--h5-path <h5_path>	moseq2-pca-apply-pca	command line option
	moseq2-pca-compute-changepoints		moseq2-pca-compute-changepoints	command line option
	moseq2-pca-train-pca	--input-dir <input_dir>	moseq2-pca-apply-pca	command line option
cluster-type <cluster_type>	moseq2-pca-apply-pca		moseq2-pca-compute-changepoints	command line option
	moseq2-pca-compute-changepoints		moseq2-pca-train-pca	command line option
	moseq2-pca-train-pca	--klags <klags>	moseq2-pca-compute-changepoints	command line option
--config-file <config_file>	moseq2-pca-apply-pca	--local-processes <local_processes>	moseq2-pca-train-pca	command line option
	moseq2-pca-compute-changepoints	--mask-height-threshold <mask_height_threshold>	moseq2-pca-train-pca	command line option
	moseq2-pca-train-pca	--mask-threshold <mask_threshold>	moseq2-pca-train-pca	command line option
--cores <cores>	moseq2-pca-apply-pca	--max-height <max_height>	moseq2-pca-train-pca	command line option
	moseq2-pca-compute-changepoints	--medfilter-space <medfilter_space>	moseq2-pca-train-pca	command line option
	moseq2-pca-train-pca	--medfilter-time <medfilter_time>	moseq2-pca-train-pca	command line option
		--memory <memory>	moseq2-pca-apply-pca	command line option
mask-cache-path <dask_cache_path>	moseq2-pca-apply-pca		moseq2-pca-compute-changepoints	command line option
	moseq2-pca-compute-changepoints		moseq2-pca-train-pca	command line option
	moseq2-pca-train-pca	--min-height <min_height>	moseq2-pca-train-pca	command line option
--detrend-window <detrend_window>	moseq2-pca-apply-pca	--missing-data	moseq2-pca-train-pca	command line option
--dims <dims>	moseq2-pca-compute-changepoints	--missing-data-iters <missing_data_iters>	moseq2-pca-train-pca	command line option
--fill-gaps <fill_gaps>	moseq2-pca-apply-pca	--neighbors <neighbors>	moseq2-pca-compute-changepoints	command line option
--fps <fps>	moseq2-pca-apply-pca	--nworkers <nworkers>	moseq2-pca-apply-pca	command line option
	moseq2-pca-compute-changepoints		moseq2-pca-compute-changepoints	command line option
--from-end	moseq2-pca-clip-scores		moseq2-pca-train-pca	command line option
--gaussfilter-space <gaussfilter_space>	moseq2-pca-train-pca			
--gaussfilter-time <gaussfilter_time>	moseq2-pca-train-pca			

--output-dir <output_dir>	moseq2-pca-apply-pca command line option	-use-fft	moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option	-visualize-results <visualize_results>	moseq2-pca-compute-changepoints command line option
	moseq2-pca-train-pca command line option		moseq2-pca-train-pca command line option
--output-file <output_file>	moseq2-pca-apply-pca command line option	--wall-time <wall_time>	moseq2-pca-apply-pca command line option
	moseq2-pca-compute-changepoints command line option		moseq2-pca-compute-changepoints command line option
	moseq2-pca-train-pca command line option		moseq2-pca-train-pca command line option
--pca-file <pca_file>	moseq2-pca-apply-pca command line option	-c	moseq2-pca-apply-pca command line option
-e-components <pca_file_components>	moseq2-pca-compute-changepoints command line option		moseq2-pca-compute-changepoints command line option
-a-file-scores <pca_file_scores>	moseq2-pca-compute-changepoints command line option		moseq2-pca-train-pca command line option
--pca-path <pca_path>	moseq2-pca-apply-pca command line option	-d	moseq2-pca-apply-pca command line option
	moseq2-pca-compute-changepoints command line option		moseq2-pca-compute-changepoints command line option
--processes <processes>	moseq2-pca-apply-pca command line option		moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option	-i	moseq2-pca-apply-pca command line option
	moseq2-pca-train-pca command line option		moseq2-pca-compute-changepoints command line option
--queue <queue>	moseq2-pca-apply-pca command line option		moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option	-k	moseq2-pca-compute-changepoints command line option
	moseq2-pca-train-pca command line option	-m	moseq2-pca-apply-pca command line option
--rank <rank>	moseq2-pca-train-pca command line option		moseq2-pca-compute-changepoints command line option
--recon-pcs <recon_pcs>	moseq2-pca-train-pca command line option		moseq2-pca-train-pca command line option
--sigma <sigma>	moseq2-pca-compute-changepoints command line option	-n	moseq2-pca-apply-pca command line option
--tailfilter-shape <tailfilter_shape>	moseq2-pca-train-pca command line option		moseq2-pca-compute-changepoints command line option
--tailfilter-size <tailfilter_size>	moseq2-pca-train-pca command line option		moseq2-pca-train-pca command line option
--threshold <threshold>	moseq2-pca-compute-changepoints command line option	-o	moseq2-pca-apply-pca command line option
--timeout <timeout>	moseq2-pca-apply-pca command line option		moseq2-pca-compute-changepoints command line option
	moseq2-pca-compute-changepoints command line option		moseq2-pca-train-pca command line option
	moseq2-pca-train-pca command line option		

`-p` moseq2-pca-apply-pca command line option
moseq2-pca-compute-changepoints command line option
moseq2-pca-train-pca command line option

`-q` moseq2-pca-apply-pca command line option
moseq2-pca-compute-changepoints command line option
moseq2-pca-train-pca command line option

`-s` moseq2-pca-compute-changepoints command line option

`-w` moseq2-pca-apply-pca command line option
moseq2-pca-compute-changepoints command line option
moseq2-pca-train-pca command line option

A

`apply_pca_command()` (in module moseq2_pca.gui)
`apply_pca_dask()` (in module moseq2_pca.pca.util)
`apply_pca_local()` (in module moseq2_pca.pca.util)
`apply_pca_wrapper()` (in module moseq2_pca.helpers.wrappers)

C

`changepoint_dist()` (in module moseq2_pca.viz)
`clean_frames()` (in module moseq2_pca.util)

CLIP_SAMPLES

moseq2-pca-clip-scores command line option
`command_with_config()` (in module moseq2_pca.util)
`compute_changepoints_command()` (in module moseq2_pca.gui)
`compute_changepoints_wrapper()` (in module moseq2_pca.helpers.wrappers)

D

`display_components()` (in module moseq2_pca.viz)

G

`gauss_smooth()` (in module moseq2_pca.util)
`gaussian_kernel1d()` (in module moseq2_pca.util)
`get_changepoints()` (in module moseq2_pca.util)
`get_changepoints_dask()` (in module moseq2_pca.pca.util)

`get_metadata_path()` (in module moseq2_pca.util)
`get_pca_yaml_data()` (in module moseq2_pca.helpers.data)
`get_rps()` (in module moseq2_pca.util)
`get_timestamp_path()` (in module moseq2_pca.util)

I

`initialize_dask()` (in module moseq2_pca.util)
`insert_nans()` (in module moseq2_pca.util)

L

`load_pcs_for_cp()` (in module moseq2_pca.helpers.data)

M

`mask_data()` (in module moseq2_pca.pca.util)

moseq2-pca-apply-pca command line option

`--chunk-size <chunk_size>`
`--cluster-type <cluster_type>`
`--config-file <config_file>`
`--cores <cores>`
`--dask-cache-path <dask_cache_path>`
`--detrend-window <detrend_window>`
`--fill-gaps <fill_gaps>`
`--fps <fps>`
`--h5-mask-path <h5_mask_path>`
`--h5-path <h5_path>`
`--input-dir <input_dir>`
`--memory <memory>`
`--nworkers <nworkers>`
`--output-dir <output_dir>`
`--output-file <output_file>`
`--pca-file <pca_file>`
`--pca-path <pca_path>`
`--processes <processes>`
`--queue <queue>`
`--timeout <timeout>`
`--wall-time <wall_time>`

`-c`
`-d`
`-i`
`-m`
`-n`
`-o`

-p
-q
-w

moseq2-pca-clip-scores command line option

--from-end
CLIP_SAMPLES
PCA_FILE

moseq2-pca-compute-changepoints command line option

--chunk-size <chunk_size>
--cluster-type <cluster_type>
--config-file <config_file>
--cores <cores>
--dask-cache-path <dask_cache_path>
--dims <dims>
--fps <fps>
--h5-mask-path <h5_mask_path>
--h5-path <h5_path>
--input-dir <input_dir>
--klags <klags>
--memory <memory>
--neighbors <neighbors>
--nworkers <nworkers>
--output-dir <output_dir>
--output-file <output_file>
--pca-file-components <pca_file_components>
--pca-file-scores <pca_file_scores>
--pca-path <pca_path>
--processes <processes>
--queue <queue>
--sigma <sigma>
--threshold <threshold>
--timeout <timeout>
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-k
-m
-n
-o
-p

-q
-s
-w

moseq2-pca-train-pca command line option

--chunk-size <chunk_size>
--cluster-type <cluster_type>
--config-file <config_file>
--cores <cores>
--dask-cache-path <dask_cache_path>
--gaussfilter-space <gaussfilter_space>
--gaussfilter-time <gaussfilter_time>
--input-dir <input_dir>
--local-processes <local_processes>
--mask-height-threshold <mask_height_threshold>
--mask-threshold <mask_threshold>
--max-height <max_height>
--medfilter-space <medfilter_space>
--medfilter-time <medfilter_time>
--memory <memory>
--min-height <min_height>
--missing-data
--missing-data-iters <missing_data_iters>
--nworkers <nworkers>
--output-dir <output_dir>
--output-file <output_file>
--processes <processes>
--queue <queue>
--rank <rank>
--recon-pcs <recon_pcs>
--tailfilter-shape <tailfilter_shape>
--tailfilter-size <tailfilter_size>
--timeout <timeout>
--use-fft
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-m
-n
-o
-p

-q

-w

moseq2_pca.gui (module)

moseq2_pca.helpers.data (module)

moseq2_pca.helpers.wrappers (module)

moseq2_pca.pca.util (module)

moseq2_pca.util (module)

moseq2_pca.viz (module)

P

PCA_FILE

moseq2-pca-clip-scores command line option

R

read_yaml() (in module moseq2_pca.util)

recursive_find_h5s() (in module moseq2_pca.util)

recursively_load_dict_contents_from_group() (in
module moseq2_pca.util)

S

scree_plot() (in module moseq2_pca.viz)

select_strel() (in module moseq2_pca.util)

setup_cp_command() (in module
moseq2_pca.helpers.data)

shutdown_dask() (in module moseq2_pca.util)

T

train_pca_command() (in module moseq2_pca.gui)

train_pca_dask() (in module moseq2_pca.pca.util)

train_pca_wrapper() (in module
moseq2_pca.helpers.wrappers)

Python Module Index

m

[moseq2_pca](#)

[moseq2_pca.gui](#)

[moseq2_pca.helpers.data](#)

[moseq2_pca.helpers.wrappers](#)

[moseq2_pca.pca.util](#)

[moseq2_pca.util](#)

[moseq2_pca.viz](#)