

Python Documentation

version

May 13, 2020

Contents

Welcome to moseq2-pca's documentation!	1
moseq2-pca	1
moseq2_pca package	1
Subpackages	1
moseq2_pca.helpers package	1
Helpers - Data Module	1
Helpers - Wrapper Module	1
moseq2_pca.pca package	2
PCA - Utilities Module	2
CLI Module	5
cli	5
apply-pca	5
clip-scores	6
compute-changepoints	6
train-pca	7
GUI Module	9
Utilities Module	9
Visualization Module	12
Indices and tables	13
Index	15
Python Module Index	21

Welcome to moseq2-pca's documentation!

moseq2-pca

moseq2_pca package

Subpackages

moseq2_pca.helpers package

Helpers - Data Module

`moseq2_pca.helpers.data.get_pca_yaml_data(pca_yaml)`

Reads PCA yaml file and returns metadata

Parameters: `pca_yaml (str)` (path to `pca.yaml`)

Returns: `use_fft (bool)` (indicates whether to use FFT) `clean_params (dict)` (dict of image filtering parameters) `mask_params (dict)` (dict of mask parameters)) `missing_data (bool)` (indicates whether to use `mask_params`)

`moseq2_pca.helpers.data.load_pcs_for_cp(pca_file_components, config_data)`

Load computed Principal Components for Model-free Changepoint Analysis.

Parameters:

- `pca_file_components (str)` (path to `pca h5` file to read PCs)
- `config_data (dict)` (config parameters)

Returns: `pca_components (str)` (path to `pca` components) `changepoint_params (dict)` (dict of relevant changepoint parameters) `cluster (dask Cluster)` (Dask Cluster object.) `client (dask Client)` (Dask Client Object) `missing_data (bool)` (Indicates whether to use `mask_params`) `mask_params (dict)` (Mask parameters to use when computing CPs)

`moseq2_pca.helpers.data.setup_cp_command(input_dir, config_data, output_dir, output_file, output_directory=None)`

Helper function for `changepoints_wrapper` to perform data-path existence checks.

Parameters:

- `input_dir (int)` (path to directory containing all `h5+yaml` files)
- `config_data (dict)` (dict of relevant PCA parameters (image filtering etc.))
- `output_dir (str)` (path to directory to store PCA data)
- `output_file (str)` (`pca` model filename)
- `output_directory (str)` (alternative `output_dir`)

Returns: `config_data (dict)` (updated `config_data` dict with the proper paths) `pca_file_components (str)` (path to trained `pca` file) `pca_file_scores (str)` (path to `pca_scores` file) `h5s (list)` (list of relevant `pca h5` files) `yamls (list)` (list of relevant `pca` metadata `yaml` files) `save_file (str)` (path to save changepoints)

Helpers - Wrapper Module

`moseq2_pca.helpers.wrappers.apply_pca_wrapper(input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)`

Wrapper function to obtain PCA Scores.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.compute_changepoints_wrapper` (input_dir, config_data, output_dir, output_file, gui=False, output_directory=None)

Wrapper function to compute model-free (PCA based) Changepoints.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.train_pca_wrapper` (input_dir, config_data, output_dir, output_file, output_directory=None, gui=False)

Wrapper function to train PCA.

Parameters:

- **input_dir (int)** (*path to directory containing all h5+yaml files*)
- **config_data (dict)** (*dict of relevant PCA parameters (image filtering etc.)*)
- **output_dir (str)** (*path to directory to store PCA data*)
- **output_file (str)** (*pca model filename*)
- **output_directory (str)** (*alternative output_dir*)
- **gui (bool)** (*indicate GUI is running*)

Returns: **config_data (dict)**

Return type: updated config_data variable to write back in GUI API

moseq2_pca.pca package

PCA - Utilities Module

`moseq2_pca.pca.util.apply_pca_dask` (pca_components, h5s, yamls, use_fft, clean_params, save_file, chunk_size, mask_params, missing_data, client, fps=30, gui=False)

"Apply" trained PCA on input frame data to obtain PCA Scores using Distributed Dask cluster.

Parameters:

- **pca_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use_fft (bool)** (indicate whether to use 2D-FFT)
- **clean_params (dict)** (dictionary containing filtering options)
- **save_file (str)** (path to pca_scores filename to save)
- **chunk_size (int)** (size of chunks to process)
- **mask_params (dict)** (dictionary of masking parameters (if missing data))
- **missing_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)

Returns:

Return type: None

`moseq2_pca.pca.util.apply_pca_local` (pca_components, h5s, yamls, use_fft, clean_params, save_file, chunk_size, mask_params, missing_data, fps=30)

“Apply” trained PCA on input frame data to obtain PCA Scores using local cluster/platform.

Parameters:

- **pca_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use_fft (bool)** (indicate whether to use 2D-FFT)
- **clean_params (dict)** (dictionary containing filtering options)
- **save_file (str)** (path to pca_scores filename to save)
- **chunk_size (int)** (size of chunks to process)
- **mask_params (dict)** (dictionary of masking parameters (if missing data))
- **missing_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)

Returns:

Return type: None

`moseq2_pca.pca.util.get_changepoints_dask` (changepoint_params, pca_components, h5s, yamls, save_file, chunk_size, mask_params, missing_data, client, fps=30, pca_scores=None, progress_bar=False, gui=False)

Computes model-free changepoints using PCs and PC Scores on distributed dask cluster.

Parameters:

- **changepoint_params (dict)** (*dict of changepoint parameters*)
- **pca_components (np.array)** (*computed principal components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **save_file (str)** (*path to save changepoint files*)
- **chunk_size (int)** (*size of chunks to process in dask.*)
- **mask_params (dict)** (*dict of missing_data mask parameters.*)
- **missing_data (bool)** (*indicate whether to use mask_params*)
- **client (dask Client)** (*initialized Dask Client object*)
- **fps (int)** (*frames per second*)
- **pca_scores (np.array)** (*computed principal component scores*)
- **progress_bar (bool)** (*display progress bar*)
- **gui (bool)** (*indicate GUI use*)

Returns:

Return type: None

`moseq2_pca.pca.util.mask_data` (`original_data`, `mask`, `new_data`)
Create a mask subregion given a boolean mask if missing data flag is used.

Parameters:

- **original_data (3d np.ndarray)** (*input frames*)
- **mask (3d boolean np.ndarray)** (*mask array*)
- **new_data (3d np.ndarray)** (*frames to use*)

Returns: output (3d np.ndarray)

Return type: masked data array

`moseq2_pca.pca.util.train_pca_dask` (`dask_array`, `clean_params`, `use_fft`, `rank`, `cluster_type`, `client`, `workers`, `cache`, `mask=None`, `iters=10`, `recon_pcs=10`, `min_height=10`, `max_height=100`)
Train PCA using dask arrays.

Parameters:

- **dask_array (dask array)** (*chunked frames to train PCA*)
- **clean_params (dict)** (*dictionary containing filtering parameters*)
- **use_fft (bool)** (*indicates whether to use 2d-FFT on images.*)
- **rank (int)** (*Matrix rank to use*)
- **cluster_type (str)** (*indicates which cluster to use.*)
- **client (Dask.Client)** (*client object to execute dask operations*)
- **workers (int)** (*number of dask workers*)
- **cache (str)** (*path to cache directory*)
- **mask (dask array)** (*dask array of masked data if missing_data parameter==True*)
- **iters (int)** (*number of SVD iterations*)
- **recon_pcs (int)** (*number of PCs to reconstruct. (if missing_data = True)*)
- **min_height (int)** (*minimum mouse height from floor in (mm)*)
- **max_height (int)** (*maximum mouse height from floor in (mm)*)

Returns: output_dict (dict)

Return type: dictionary containing PCA training results.

CLI Module

cli

```
cli [OPTIONS] COMMAND [ARGS]...
```

apply-pca

```
cli apply-pca [OPTIONS]
```

Options

-i, --input-dir <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

--cluster-type <cluster_type>
Cluster type [default: local]
Options: local|slurm|nodask

-o, --output-dir <output_dir>
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca]

--output-file <output_file>
Name of h5 file for storing pca results [default: pca_scores]

--h5-path <h5_path>
Path to data in h5 files [default: /frames]

--h5-mask-path <h5_mask_path>
Path to log-likelihood mask in h5 files [default: /frames_mask]

--pca-path <pca_path>
Path to pca components [default: /components]

--pca-file <pca_file>
Path to PCA results

--chunk-size <chunk_size>
Number of frames per chunk [default: 4000]

--fill-gaps <fill_gaps>
Fill dropped frames with nans [default: True]

--fps <fps>
Fps (only used if no timestamps found) [default: 30]

--detrend-window <detrend_window>
Length of detrend window (in seconds, 0 for no detrending) [default: 0]

--config-file <config_file>
Path to configuration file

-d, --dask-cache-path <dask_cache_path>
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

-q, --queue <queue>
Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
Number of workers [default: 10]

-c, --cores <cores>
Number of cores per worker [default: 1]

-p, --processes <processes>
Number of processes to run on each worker [default: 1]

-m, --memory <memory>

Welcome to moseq2-pca's documentation!

RAM usage per workers [default: 15GB]

-w, --wall-time <wall_time>

Wall time for workers [default: 06:00:00]

--timeout <timeout>

Time to wait for workers to initialize before proceeding (minutes) [default: 5]

clip-scores

Clips PCA scores from the beginning or end

Args:

pca_file (string): Path to PCA scores **clip_samples** (int): number of samples to clip from beginning or end
from_end (bool): if true clip from end rather than beginning

Note that scores are modified *in place*.

```
cli clip-scores [OPTIONS] PCA_FILE CLIP_SAMPLES
```

Options

--from-end

[default: False]

Arguments

PCA_FILE

Required argument

CLIP_SAMPLES

Required argument

compute-changepoints

```
cli compute-changepoints [OPTIONS]
```

Options

-i, --input-dir <input_dir>

Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

-o, --output-dir <output_dir>

Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca/]

--output-file <output_file>

Name of h5 file for storing pca results [default: changepoints]

--cluster-type <cluster_type>

Cluster type [default: local]

Options: local|slurm

--pca-file-components <pca_file_components>

Path to PCA components

--pca-file-scores <pca_file_scores>

Path to PCA results

--pca-path <pca_path>

Path to pca components [default: /components]

--neighbors <neighbors>

Neighbors to use for peak identification [default: 1]

--threshold <threshold>

Peak threshold to use for changepoints [default: 0.5]

-k, --klags <klags>

Lag to use for derivative calculation [default: 6]

Welcome to moseq2-pca's documentation!

-s, --sigma <sigma>
Standard deviation of gaussian smoothing filter [default: 3.5]

-d, --dims <dims>
Number of random projections to use [default: 300]

--fps <fps>
Fps (only used if no timestamps found) [default: 30]

--h5-path <h5_path>
Path to data in h5 files [default: /frames]

--h5-mask-path <h5_mask_path>
Path to log-likelihood mask in h5 files [default: /frames_mask]

--chunk-size <chunk_size>
Number of frames per chunk [default: 4000]

--config-file <config_file>
Path to configuration file

--dask-cache-path <dask_cache_path>
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

--visualize-results <visualize_results>
Visualize results [default: True]

-q, --queue <queue>
Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
Number of workers [default: 10]

-c, --cores <cores>
Number of cores per worker [default: 1]

-p, --processes <processes>
Number of processes to run on each worker [default: 1]

-m, --memory <memory>
RAM usage per workers [default: 15GB]

-w, --wall-time <wall_time>
Wall time for workers [default: 06:00:00]

--timeout <timeout>
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

train-pca

```
cli train-pca [OPTIONS]
```

Options

-i, --input-dir <input_dir>
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

--cluster-type <cluster_type>
Cluster type [default: local]
Options: local|slurm

-o, --output-dir <output_dir>
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/_pca]

--gaussfilter-space <gaussfilter_space>
Spatial filter for data (Gaussian) [default: 1.5, 1]

--gaussfilter-time <gaussfilter_time>
Temporal filter for data (Gaussian) [default: 0]

Welcome to moseq2-pca's documentation!

```
--medfilter-space <medfilter_space>
  Median spatial filter [default: 0]

--medfilter-time <medfilter_time>
  Median temporal filter [default: 0]

--missing-data
  Use missing data PCA [default: False]

--missing-data-iters <missing_data_iters>
  Missing data PCA iterations [default: 10]

--mask-threshold <mask_threshold>
  Threshold for mask (missing data only) [default: -16]

--mask-height-threshold <mask_height_threshold>
  Threshold for mask based on floor height [default: 5]

--min-height <min_height>
  Min mouse height from floor (mm) [default: 10]

--max-height <max_height>
  Max mouse height from floor (mm) [default: 100]

--tailfilter-size <tailfilter_size>
  Tail filter size [default: 9, 9]

--tailfilter-shape <tailfilter_shape>
  Tail filter shape [default: ellipse]

--use-fft
  Use 2D fft [default: False]

--recon-pcs <recon_pcs>
  Number of PCs to use for missing data reconstruction [default: 10]

--rank <rank>
  Rank for compressed SVD (generally>>nPCS) [default: 50]

--output-file <output_file>
  Name of h5 file for storing pca results [default: pca]

--chunk-size <chunk_size>
  Number of frames per chunk [default: 4000]

--visualize-results <visualize_results>
  Visualize results [default: True]

--config-file <config_file>
  Path to configuration file

-d, --dask-cache-path <dask_cache_path>
  Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

--local-processes <local_processes>
  Use processes with local scheduler [default: True]

-q, --queue <queue>
  Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
  Number of workers [default: 10]

-c, --cores <cores>
  Number of cores per worker [default: 1]

-p, --processes <processes>
  Number of processes to run on each worker [default: 1]

-m, --memory <memory>
  Total RAM usage per worker [default: 15GB]

-w, --wall-time <wall_time>
```

Wall time for workers [default: 06:00:00]

--timeout <timeout>

Time to wait for workers to initialize before proceeding (minutes) [default: 5]

GUI Module

`moseq2_pca.gui.apply_pca_command` (input_dir, index_file, config_file, output_dir, output_file, output_directory=None)

Compute PCA Scores given trained PCA using Jupyter Notebook.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **index_file (str)** (path to index file.)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.compute_changepoints_command` (input_dir, config_file, output_dir, output_file, output_directory=None)

Compute Changepoint distribution using Jupyter Notebook.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns: (str)

Return type: success string.

`moseq2_pca.gui.train_pca_command` (input_dir, config_file, output_dir, output_file, output_directory=None)

Train PCA through Jupyter notebook, and updates config file.

Parameters:

- **input_dir (str)** (path to directory containing training data)
- **config_file (str)** (path to config file)
- **output_dir (str)** (path to output pca directory)
- **output_file (str)** (name of output pca file.)
- **output_directory (str)** (alternative output directory path)

Returns:

Return type: None

Utilities Module

`moseq2_pca.util.clean_frames` (frames, medfilter_space=None, gaussfilter_space=None, medfilter_time=None, gaussfilter_time=None, detrend_time=None, tailfilter=None, tail_threshold=5)

Filters spatial/temporal noise from frames using Median and Gaussian filters, given kernel sizes for each respective requested filter.

Parameters:

- **frames (3D numpy array)** (*frames to filter.*)
- **medfilter_space (list)** (*median spatial filter kernel.*)
- **gaussfilter_space (list)** (*gaussian spatial filter kernel.*)
- **medfilter_time (list)** (*median temporal filter.*)
- **gaussfilter_time (list)** (*gaussian temporal filter.*)
- **detrend_time (int)** (*number of frames to lag for.*)
- **tailfilter (int)** (*size of tail-filter kernel.*)
- **tail_threshold (int)** (*threshold value to use for tail filtering*)

Returns: **out (3D numpy array)**

Return type: filtered frames.

`moseq2_pca.util.command_with_config` (config_file_param_name)

`moseq2_pca.util.gauss_smooth` (signal, win_length=None, sig=1.5, kernel=None)
Perform Gaussian Smoothing on a 1D signal.

Parameters:

- **signal (1d numpy array)** (*signal to perform smoothing*)
- **win_length (int)** (*window_size for gaussian kernel filter*)
- **sig (float)** (*variance of 1d gaussian kernel.*)
- **kernel (tuple)** (*kernel size to use for smoothing*)

Returns: **result (1d numpy array)**

Return type: smoothed signal

`moseq2_pca.util.gaussian_kernel1d` (n=None, sig=3)
Get 1D gaussian kernel.

Parameters:

- **n (int)** (*number of points to use.*)
- **sig (int)** (*variance of kernel to use.*)

Returns: **kernel (1d array)**

Return type: 1D numpy kernel.

`moseq2_pca.util.get_changepoints` (scores, k=5, sigma=3, peak_height=0.5, peak_neighbors=1, baseline=True, timestamps=None)
Compute changepoints distribution and CP Curve.

Parameters:

- **scores (3D numpy array)** (*nframes * r * c*)
- **k (int)** (*klags - Lag to use for derivative calculation.*)
- **sigma (int)** (*Standard deviation of gaussian smoothing filter.*)
- **peak_height (float)** (*user-defined peak Changepoint length.*)
- **peak_neighbors (int)** (*number of peaks in the CP curve.*)
- **baseline (bool)** (*normalize data.*)
- **timestamps (array)** (*loaded timestamps.*)

Returns: **cps (2D numpy array)** (*array of values for CP curve*) **normed_df (1D numpy array)** (*array of values for bar plot*)

`moseq2_pca.util.get_metadata_path` (h5file)
Return path within h5 file that contains the kinect extraction metadata.

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to acquisition metadata within h5 file.

Welcome to moseq2-pca's documentation!

```
moseq2_pca.util.get_rps(frames, rps=600, normalize=True)
```

Get random projections of frames.

Parameters:

- **frames (2D or 3D numpy array)** (*Frames to get dimensions from.*)
- **rps (int)** (*Number of random projections.*)
- **normalize (bool)** (*indicates whether to normalize frames.*)

Returns: **rproj (2D or 3D numpy array)**

Return type: Computed random projections with same shape as frames

```
moseq2_pca.util.get_timestamp_path(h5file)
```

Return path within h5 file that contains the kinect timestamps

Parameters: **h5file (str)** (*path to h5 file.*)

Returns: **(str)**

Return type: path to metadata timestamps within h5 file

```
moseq2_pca.util.initialize_dask(nworkers=50, processes=1, memory='4GB', cores=1,
wall_time='01:00:00', queue='debug', local_processes=False, cluster_type='local',
scheduler='distributed', timeout=10, cache_path='/Users/aymanzeine/moseq2_pca', **kwargs)
```

Initialize dask client, cluster, workers, etc.

Parameters:

- **nworkers (int)** (*number of dask workers to initialize*)
- **processes (int)** (*number of processes per worker*)
- **memory (str)** (*amount of memory to allocate to dask cluster*)
- **cores (int)** (*number of cores to use.*)
- **wall_time (str)** (*amount of time to allow program to run*)
- **queue (str)** (*logging mode*)
- **local_processes (bool)** (*indicate whether the processes are local*)
- **cluster_type (str)** (*indicate what cluster to use*)
- **scheduler (str)** (*indicate what scheduler to use*)
- **timeout (int)** (*number of worker timeouts to allow*)
- **cache_path (str or Pathlike)** (*path to store cached data*)
- **kwargs** (*extra keyword arguments*)

Returns: **client (dask Client)** (*initialized Client*) **cluster (dask Cluster)** (*initialized Cluster*) **workers (dask Workers)** (*initialized workers or None if cluster_type = 'local'*) **cache (dask Chest)** (*initialized Chest (cache) object pointing to given cache path*)

```
moseq2_pca.util.insert_nans(timestamps, data, fps=30)
```

Fills NaN values with 0 in timestamps.

Parameters:

- **timestamps (1D array)** (*timestamp time-strs*)
- **data (1D array)** (*timestamp values*)
- **fps (int)** (*frames per second*)

Returns: **filled_data (1D array)** (*filled missing timestamp values.*) **data_idx (1D array)** (*indices of inserted 0s*) **filled_timestamps (1D array)** (*filled timestamp-strs*)

```
moseq2_pca.util.read_yaml(yaml_file)
```

Reads yaml file and returns dictionary representation of file contents.

Parameters: **yaml_file (str)** (*path to yaml file*)

Returns: **return_dict (dict)**

Return type: dict of yaml file contents

```
moseq2_pca.util.recursive_find_h5s  
(root_dir='/Users/aymanzeine/Desktop/moseq/moseq2-pca/docs', ext='.h5',  
yaml_string='{}.yaml')  
Recursively find h5 files, along with yaml files with the same basename
```

Parameters:

- **root_dir (str or os.Pathlike)** (path to directory to start recursive search)
- **ext (str)** (extension to search for, e.g. .h5)
- **yaml_string (str)** (a format to use to name yaml files)

Returns: **h5s (list)** (list of h5 file paths) **dicts (list)** (list of dicts containing metadata file contents)
yamls (list) (list of yaml file paths)

```
moseq2_pca.util.recursively_load_dict_contents_from_group (h5file, path)  
Reads all contents from h5 and returns them in a nested dict object.
```

Parameters:

- **h5file (str)** (path to h5 file)
- **path (str)** (path to group within h5 file)

Returns: **ans (dict)**

Return type: dictionary of all h5 group contents

```
moseq2_pca.util.select_strel (string='e', size=(10,10))  
Selects Structuring Element Shape
```

Parameters:

- **string (str)** (e for Ellipse, r for Rectangle)
- **size (tuple)** (size of StructuringElement)

Returns: **strel (cv2.StructuringElement)**

Return type: returned StructuringElement with specified size.

```
moseq2_pca.util.shutdown_dask (scheduler)  
Graceful shutdown dask scheduler.  
https://github.com/dask/distributed/issues/1703#issuecomment-361291492
```

Parameters: **scheduler (dask Scheduler)** (scheduler to shutdown.)

Returns:

Return type: None

Visualization Module

```
moseq2_pca.viz.changepoint_dist (cps, headless=False)  
Creates bar plot describing computed Changepoint Distribution.
```

Parameters:

- **cps (np.ndarray)** (changepoints to graph)
- **headless (bool)** (trim first element in PC list)

Returns: **plt (plt.figure)** (figure to save/graph) **ax (plt.ax)** (figure axis variable)

```
moseq2_pca.viz.display_components (components, cmap='gray', headless=False)  
Creates grid of computed Principal Components.
```

Parameters:

- **components (2D np.ndarray)** (components to graph)
- **cmap (str)** (color map to use)
- **headless (bool)** (trim first element in PC list)

Returns: **plt (plt.figure)** (figure to save/graph) **ax (plt.ax)** (figure axis variable)

```
moseq2_pca.viz.scree_plot (explained_variance_ratio, headless=False)  
Creates Scree plot describing principal components.
```


Parameters:

- **explained_variance_ratio** (1D `np.array`) (*explained variance ratio of each principal component*)
- **headless** (`bool`) (*trim first element in PC list*)

Returns: `plt` (`plt.figure`)

Return type: figure to save/graph

Indices and tables

- `genindex`
- `modindex`
- `search`

Index

Symbols

			--gaussfilter-time <gaussfilter_time>	cli-train-pca command line option
			--h5-mask-path <h5_mask_path>	cli-apply-pca line option
--chunk-size <chunk_size>	cli-apply-pca line option	command		cli-compute-changepoints command line option
	cli-compute-changepoints command line option		--h5-path <h5_path>	cli-apply-pca line option
	cli-train-pca line option	command		cli-compute-changepoints command line option
--cluster-type <cluster_type>	cli-apply-pca line option	command	--input-dir <input_dir>	cli-apply-pca line option
	cli-compute-changepoints command line option			cli-compute-changepoints command line option
	cli-train-pca line option	command		cli-train-pca line option
--config-file <config_file>	cli-apply-pca line option	command	--klags <klags>	cli-compute-changepoints command line option
	cli-compute-changepoints command line option		local-processes <local_processes>	cli-train-pca command line option
	cli-train-pca line option	command	--mask-height-threshold <mask_height_threshold>	cli-train-pca command line option
--cores <cores>	cli-apply-pca line option	command		
	cli-compute-changepoints command line option		--mask-threshold <mask_threshold>	cli-train-pca command line option
	cli-train-pca line option	command	--max-height <max_height>	cli-train-pca command line option
--dask-cache-path <dask_cache_path>	cli-apply-pca line option	command		
	cli-compute-changepoints command line option		--medfilter-space <medfilter_space>	cli-train-pca command line option
	cli-train-pca line option	command	--medfilter-time <medfilter_time>	cli-train-pca command line option
--detrend-window <detrend_window>	cli-apply-pca command line option		--memory <memory>	cli-apply-pca line option
	cli-compute-changepoints command line option			cli-compute-changepoints command line option
--dims <dims>	cli-apply-pca command line option			cli-train-pca command line option
--fill-gaps <fill_gaps>	cli-apply-pca command line option		--min-height <min_height>	cli-train-pca command line option
--fps <fps>	cli-apply-pca line option	command		
	cli-compute-changepoints command line option		--missing-data	cli-train-pca command line option
--from-end	cli-clip-scores line option	command	--missing-data-iters <missing_data_iters>	cli-train-pca command line option
--gaussfilter-space <gaussfilter_space>	cli-train-pca command line option		--neighbors <neighbors>	cli-compute-changepoints command line option

--nworkers <nworkers>	cli-apply-pca line option	command	--threshold <threshold>	cli-compute-changepoints command line option	
	cli-compute-changepoints command line option		--timeout <timeout>	cli-apply-pca line option	command
	cli-train-pca line option	command		cli-compute-changepoints command line option	
--output-dir <output_dir>	cli-apply-pca line option	command		cli-train-pca line option	command
	cli-compute-changepoints command line option		--use-fft	cli-train-pca line option	command
	cli-train-pca line option	command	--visualize-results <visualize_results>	cli-compute-changepoints command line option	
--output-file <output_file>	cli-apply-pca line option	command		cli-train-pca line option	command
	cli-compute-changepoints command line option		--wall-time <wall_time>	cli-apply-pca line option	command
	cli-train-pca line option	command		cli-compute-changepoints command line option	
--pca-file <pca_file>	cli-apply-pca command line option			cli-train-pca line option	command
-file-components <pca_file_components>	cli-compute-changepoints command line option			cli-apply-pca line option	command
-pca-file-scores <pca_file_scores>	cli-compute-changepoints command line option			cli-compute-changepoints command line option	
--pca-path <pca_path>	cli-apply-pca line option	command		cli-train-pca line option	command
	cli-compute-changepoints command line option		-d	cli-apply-pca line option	command
--processes <processes>	cli-apply-pca line option	command		cli-compute-changepoints command line option	
	cli-compute-changepoints command line option			cli-train-pca line option	command
	cli-train-pca line option	command	-i	cli-apply-pca line option	command
--queue <queue>	cli-apply-pca line option	command		cli-compute-changepoints command line option	
	cli-compute-changepoints command line option			cli-train-pca line option	command
	cli-train-pca line option	command	-k	cli-compute-changepoints command line option	
--rank <rank>	cli-train-pca line option	command	-m	cli-apply-pca line option	command
--recon-pcs <recon_pcs>	cli-train-pca command line option			cli-compute-changepoints command line option	
--sigma <sigma>	cli-compute-changepoints command line option			cli-train-pca line option	command
--tailfilter-shape <tailfilter_shape>	cli-train-pca command line option		-n	cli-apply-pca line option	command
	cli-train-pca command line option			cli-compute-changepoints command line option	
--tailfilter-size <tailfilter_size>	cli-train-pca command line option			cli-train-pca line option	command

-o	cli-apply-pca	command line option	--input-dir <input_dir>
	cli-compute-changepoints	command line option	--memory <memory>
	cli-train-pca	command line option	--nworkers <nworkers>
-p	cli-apply-pca	command line option	--output-dir <output_dir>
	cli-compute-changepoints	command line option	--output-file <output_file>
	cli-train-pca	command line option	--pca-file <pca_file>
-q	cli-apply-pca	command line option	--pca-path <pca_path>
	cli-compute-changepoints	command line option	--processes <processes>
	cli-train-pca	command line option	--queue <queue>
-s	cli-apply-pca	command line option	--timeout <timeout>
	cli-compute-changepoints	command line option	--wall-time <wall_time>
	cli-train-pca	command line option	-c
-w	cli-apply-pca	command line option	-d
	cli-compute-changepoints	command line option	-i
	cli-train-pca	command line option	-m
			-n
			-o
			-p
			-q
			-w

A

[apply_pca_command\(\)](#) (in module `moseq2_pca.gui`)
[apply_pca_dask\(\)](#) (in module `moseq2_pca.pca.util`)
[apply_pca_local\(\)](#) (in module `moseq2_pca.pca.util`)
[apply_pca_wrapper\(\)](#) (in module `moseq2_pca.helpers.wrappers`)

C

[changepoint_dist\(\)](#) (in module `moseq2_pca.viz`)
[clean_frames\(\)](#) (in module `moseq2_pca.util`)

cli-apply-pca command line option

--chunk-size <chunk_size>
 --cluster-type <cluster_type>
 --config-file <config_file>
 --cores <cores>
 --dask-cache-path <dask_cache_path>
 --detrend-window <detrend_window>
 --fill-gaps <fill_gaps>
 --fps <fps>
 --h5-mask-path <h5_mask_path>
 --h5-path <h5_path>

cli-clip-scores command line option

--from-end
 CLIP_SAMPLES
 PCA_FILE

cli-compute-changepoints command line option

--chunk-size <chunk_size>
 --cluster-type <cluster_type>
 --config-file <config_file>
 --cores <cores>
 --dask-cache-path <dask_cache_path>
 --dims <dims>
 --fps <fps>
 --h5-mask-path <h5_mask_path>
 --h5-path <h5_path>
 --input-dir <input_dir>
 --klags <klags>
 --memory <memory>
 --neighbors <neighbors>
 --nworkers <nworkers>
 --output-dir <output_dir>
 --output-file <output_file>
 --pca-file-components <pca_file_components>

```

--pca-file-scores <pca_file_scores>
--pca-path <pca_path>
--processes <processes>
--queue <queue>
--sigma <sigma>
--threshold <threshold>
--timeout <timeout>
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-k
-m
-n
-o
-p
-q
-s
-w

```

cli-train-pca command line option

```

--chunk-size <chunk_size>
--cluster-type <cluster_type>
--config-file <config_file>
--cores <cores>
--dask-cache-path <dask_cache_path>
--gaussfilter-space <gaussfilter_space>
--gaussfilter-time <gaussfilter_time>
--input-dir <input_dir>
--local-processes <local_processes>
--mask-height-threshold <mask_height_threshold>
--mask-threshold <mask_threshold>
--max-height <max_height>
--medfilter-space <medfilter_space>
--medfilter-time <medfilter_time>
--memory <memory>
--min-height <min_height>
--missing-data
--missing-data-iters <missing_data_iters>
--nworkers <nworkers>
--output-dir <output_dir>
--output-file <output_file>

```

```

--processes <processes>
--queue <queue>
--rank <rank>
--recon-pcs <recon_pcs>
--tailfilter-shape <tailfilter_shape>
--tailfilter-size <tailfilter_size>
--timeout <timeout>
--use-fft
--visualize-results <visualize_results>
--wall-time <wall_time>
-c
-d
-i
-m
-n
-o
-p
-q
-w

```

CLIP_SAMPLES

cli-clip-scores command line option

command_with_config() (in module moseq2_pca.util)
 compute_changepoints_command() (in module moseq2_pca.gui)
 compute_changepoints_wrapper() (in module moseq2_pca.helpers.wrappers)

D

display_components() (in module moseq2_pca.viz)

G

gauss_smooth() (in module moseq2_pca.util)
 gaussian_kernel1d() (in module moseq2_pca.util)
 get_changepoints() (in module moseq2_pca.util)
 get_changepoints_dask() (in module moseq2_pca.pca.util)
 get_metadata_path() (in module moseq2_pca.util)
 get_pca_yaml_data() (in module moseq2_pca.helpers.data)
 get_rps() (in module moseq2_pca.util)
 get_timestamp_path() (in module moseq2_pca.util)

I

initialize_dask() (in module moseq2_pca.util)

[insert_nans\(\)](#) (in module [moseq2_pca.util](#))

L

[load_pcs_for_cp\(\)](#) (in module [moseq2_pca.helpers.data](#))

M

[mask_data\(\)](#) (in module [moseq2_pca.pca.util](#))

[moseq2_pca.gui](#) (module)

[moseq2_pca.helpers.data](#) (module)

[moseq2_pca.helpers.wrappers](#) (module)

[moseq2_pca.pca.util](#) (module)

[moseq2_pca.util](#) (module)

[moseq2_pca.viz](#) (module)

P

PCA_FILE

[cli-clip-scores](#) command line option

R

[read_yaml\(\)](#) (in module [moseq2_pca.util](#))

[recursive_find_h5s\(\)](#) (in module [moseq2_pca.util](#))

[recursively_load_dict_contents_from_group\(\)](#) (in module [moseq2_pca.util](#))

S

[scree_plot\(\)](#) (in module [moseq2_pca.viz](#))

[select_strel\(\)](#) (in module [moseq2_pca.util](#))

[setup_cp_command\(\)](#) (in module [moseq2_pca.helpers.data](#))

[shutdown_dask\(\)](#) (in module [moseq2_pca.util](#))

T

[train_pca_command\(\)](#) (in module [moseq2_pca.gui](#))

[train_pca_dask\(\)](#) (in module [moseq2_pca.pca.util](#))

[train_pca_wrapper\(\)](#) (in module [moseq2_pca.helpers.wrappers](#))

Python Module Index

m

[moseq2_pca](#)

[moseq2_pca.gui](#)

[moseq2_pca.helpers.data](#)

[moseq2_pca.helpers.wrappers](#)

[moseq2_pca.pca.util](#)

[moseq2_pca.util](#)

[moseq2_pca.viz](#)