

# Python Documentation

version

June 12, 2020



# Contents

<b>Welcome to moseq2-pca's documentation!</b>	<b>1</b>
moseq2_pca package	1
CLI Module	1
moseq2-pca	1
apply-pca	1
clip-scores	2
compute-changepoints	2
train-pca	3
version	5
GUI Module	5
Utilities Module	6
Visualization Module	8
Subpackages	9
moseq2_pca.helpers package	9
Helpers - Data Module	9
Helpers - Wrapper Module	10
moseq2_pca.pca package	11
PCA - Utilities Module	11
<b>Indices and tables</b>	<b>14</b>
<b>Index</b>	<b>15</b>
<b>Python Module Index</b>	<b>21</b>



# Welcome to moseq2-pca's documentation!

## moseq2\_pca package

### CLI Module

#### *moseq2-pca*

```
moseq2-pca [OPTIONS] COMMAND [ARGS]...
```

#### *apply-pca*

Computes PCA Scores of extraction data given a pre-trained PCA

```
moseq2-pca apply-pca [OPTIONS]
```

#### Options

- i, --input-dir** <input\_dir>  
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]
- cluster-type** <cluster\_type>  
Cluster type [default: local]  
**Options:** local|slurm|nodask
- o, --output-dir** <output\_dir>  
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/\_pca]
- output-file** <output\_file>  
Name of h5 file for storing pca results [default: pca\_scores]
- h5-path** <h5\_path>  
Path to data in h5 files [default: /frames]
- h5-mask-path** <h5\_mask\_path>  
Path to log-likelihood mask in h5 files [default: /frames\_mask]
- pca-path** <pca\_path>  
Path to pca components [default: /components]
- pca-file** <pca\_file>  
Path to PCA results
- chunk-size** <chunk\_size>  
Number of frames per chunk [default: 4000]
- fill-gaps** <fill\_gaps>  
Fill dropped frames with nans [default: True]
- fps** <fps>  
Fps (only used if no timestamps found) [default: 30]
- detrend-window** <detrend\_window>  
Length of detrend window (in seconds, 0 for no detrending) [default: 0]
- config-file** <config\_file>  
Path to configuration file
- d, --dask-cache-path** <dask\_cache\_path>  
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2\_pca]
- q, --queue** <queue>  
Cluster queue/partition for submitting jobs [default: debug]
- n, --nworkers** <nworkers>

Welcome to moseq2-pca's documentation!

Number of workers [default: 10]

**-c, --cores** <cores>  
Number of cores per worker [default: 1]

**-p, --processes** <processes>  
Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>  
RAM usage per workers [default: 15GB]

**-w, --wall-time** <wall\_time>  
Wall time for workers [default: 06:00:00]

**--timeout** <timeout>  
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

## **clip-scores**

Clips specified number of frames from PCA scores at the beginning or end

```
moseq2-pca clip-scores [OPTIONS] PCA_FILE CLIP_SAMPLES
```

### Options

**--from-end**  
[default: False]

### Arguments

**PCA\_FILE**  
Required argument

**CLIP\_SAMPLES**  
Required argument

## **compute-changepoints**

Computes the Model-Free Syllable Changepoints based on the PCA/PCA\_Scores

```
moseq2-pca compute-changepoints [OPTIONS]
```

### Options

**-i, --input-dir** <input\_dir>  
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

**-o, --output-dir** <output\_dir>  
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/\_pca/]

**--output-file** <output\_file>  
Name of h5 file for storing pca results [default: changepoints]

**--cluster-type** <cluster\_type>  
Cluster type [default: local]  
**Options:** local|slurm

**--pca-file-components** <pca\_file\_components>  
Path to PCA components

**--pca-file-scores** <pca\_file\_scores>  
Path to PCA results

**--pca-path** <pca\_path>  
Path to pca components [default: /components]

**--neighbors** <neighbors>  
Neighbors to use for peak identification [default: 1]

**--threshold** <threshold>

Welcome to moseq2-pca's documentation!

Peak threshold to use for changepoints [default: 0.5]

**-k, --klags** <klags>  
Lag to use for derivative calculation [default: 6]

**-s, --sigma** <sigma>  
Standard deviation of gaussian smoothing filter [default: 3.5]

**-d, --dims** <dims>  
Number of random projections to use [default: 300]

**--fps** <fps>  
Fps (only used if no timestamps found) [default: 30]

**--h5-path** <h5\_path>  
Path to data in h5 files [default: /frames]

**--h5-mask-path** <h5\_mask\_path>  
Path to log-likelihood mask in h5 files [default: /frames\_mask]

**--chunk-size** <chunk\_size>  
Number of frames per chunk [default: 4000]

**--config-file** <config\_file>  
Path to configuration file

**--dask-cache-path** <dask\_cache\_path>  
Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2\_pca]

**--visualize-results** <visualize\_results>  
Visualize results [default: True]

**-q, --queue** <queue>  
Cluster queue/partition for submitting jobs [default: debug]

**-n, --nworkers** <nworkers>  
Number of workers [default: 10]

**-c, --cores** <cores>  
Number of cores per worker [default: 1]

**-p, --processes** <processes>  
Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>  
RAM usage per workers [default: 15GB]

**-w, --wall-time** <wall\_time>  
Wall time for workers [default: 06:00:00]

**--timeout** <timeout>  
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

### ***train-pca***

Trains PCA on all extracted results (h5 files) in input directory

```
moseq2-pca train-pca [OPTIONS]
```

#### Options

**-i, --input-dir** <input\_dir>  
Directory to find h5 files [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs]

**--cluster-type** <cluster\_type>  
Cluster type [default: local]

**Options:** local|slurm

**-o, --output-dir** <output\_dir>  
Directory to store results [default: /Users/aymanzeine/Desktop/moseq/moseq2-pca/docs/\_pca]

Welcome to moseq2-pca's documentation!

```
--h5-path <h5_path>
  Path to data in h5 files [default: /frames]

--h5-mask-path <h5_mask_path>
  Path to log-likelihood mask in h5 files [default: /frames_mask]

--gaussfilter-space <gaussfilter_space>
  Spatial filter for data (Gaussian) [default: 1.5, 1]

--gaussfilter-time <gaussfilter_time>
  Temporal filter for data (Gaussian) [default: 0]

--medfilter-space <medfilter_space>
  Median spatial filter [default: 0]

--medfilter-time <medfilter_time>
  Median temporal filter [default: 0]

--missing-data
  Use missing data PCA [default: False]

--missing-data-iters <missing_data_iters>
  Missing data PCA iterations [default: 10]

--mask-threshold <mask_threshold>
  Threshold for mask (missing data only) [default: -16]

--mask-height-threshold <mask_height_threshold>
  Threshold for mask based on floor height [default: 5]

--min-height <min_height>
  Min mouse height from floor (mm) [default: 10]

--max-height <max_height>
  Max mouse height from floor (mm) [default: 100]

--tailfilter-size <tailfilter_size>
  Tail filter size [default: 9, 9]

--tailfilter-shape <tailfilter_shape>
  Tail filter shape [default: ellipse]

--use-fft
  Use 2D fft [default: False]

--recon-pcs <recon_pcs>
  Number of PCs to use for missing data reconstruction [default: 10]

--rank <rank>
  Rank for compressed SVD (generally>>nPCS) [default: 50]

--output-file <output_file>
  Name of h5 file for storing pca results [default: pca]

--chunk-size <chunk_size>
  Number of frames per chunk [default: 4000]

--visualize-results <visualize_results>
  Visualize results [default: True]

--config-file <config_file>
  Path to configuration file

-d, --dask-cache-path <dask_cache_path>
  Path to spill data to disk for dask local scheduler [default: /Users/aymanzeine/moseq2_pca]

--local-processes <local_processes>
  Use processes with local scheduler [default: True]

-q, --queue <queue>
  Cluster queue/partition for submitting jobs [default: debug]

-n, --nworkers <nworkers>
```



Welcome to moseq2-pca's documentation!

Number of workers [default: 10]

**-c, --cores** <cores>  
Number of cores per worker [default: 1]

**-p, --processes** <processes>  
Number of processes to run on each worker [default: 1]

**-m, --memory** <memory>  
Total RAM usage per worker [default: 15GB]

**-w, --wall-time** <wall\_time>  
Wall time for workers [default: 06:00:00]

**--timeout** <timeout>  
Time to wait for workers to initialize before proceeding (minutes) [default: 5]

## version

Print version number

```
moseq2-pca version [OPTIONS]
```

## GUI Module

`moseq2_pca.gui.apply_pca_command` (input\_dir, index\_file, config\_file, output\_dir, output\_file, output\_directory=None)

Compute PCA Scores given trained PCA using Jupyter Notebook.

**Parameters:**

- **input\_dir (str)** (path to directory containing training data)
- **index\_file (str)** (path to index file.)
- **config\_file (str)** (path to config file)
- **output\_dir (str)** (path to output pca directory)
- **output\_file (str)** (name of output pca file.)
- **output\_directory (str)** (alternative output directory path)

**Returns:** (str)

**Return type:** success string.

`moseq2_pca.gui.compute_changepoints_command` (input\_dir, config\_file, output\_dir, output\_file, output\_directory=None)

Compute Changepoint distribution using Jupyter Notebook.

**Parameters:**

- **input\_dir (str)** (path to directory containing training data)
- **config\_file (str)** (path to config file)
- **output\_dir (str)** (path to output pca directory)
- **output\_file (str)** (name of output pca file.)
- **output\_directory (str)** (alternative output directory path)

**Returns:** (str)

**Return type:** success string.

`moseq2_pca.gui.train_pca_command` (input\_dir, config\_file, output\_dir, output\_file, output\_directory=None)

Train PCA through Jupyter notebook, and updates config file.

**Parameters:**

- **input\_dir (str)** (path to directory containing training data)
- **config\_file (str)** (path to config file)
- **output\_dir (str)** (path to output pca directory)
- **output\_file (str)** (name of output pca file.)
- **output\_directory (str)** (alternative output directory path)

**Returns:**

**Return type:** None

## Utilities Module

`moseq2_pca.util.clean_frames` (frames, medfilter\_space=None, gaussfilter\_space=None, medfilter\_time=None, gaussfilter\_time=None, detrend\_time=None, tailfilter=None, tail\_threshold=5)

Filters spatial/temporal noise from frames using Median and Gaussian filters, given kernel sizes for each respective requested filter.

**Parameters:**

- **frames (3D numpy array)** (frames to filter.)
- **medfilter\_space (list)** (median spatial filter kernel.)
- **gaussfilter\_space (list)** (gaussian spatial filter kernel.)
- **medfilter\_time (list)** (median temporal filter.)
- **gaussfilter\_time (list)** (gaussian temporal filter.)
- **detrend\_time (int)** (number of frames to lag for.)
- **tailfilter (int)** (size of tail-filter kernel.)
- **tail\_threshold (int)** (threshold value to use for tail filtering)

**Returns:** out (3D numpy array)

**Return type:** filtered frames.

`moseq2_pca.util.command_with_config` (config\_file\_param\_name)

`moseq2_pca.util.gauss_smooth` (signal, win\_length=None, sig=1.5, kernel=None)  
Perform Gaussian Smoothing on a 1D signal.

**Parameters:**

- **signal (1d numpy array)** (signal to perform smoothing)
- **win\_length (int)** (window\_size for gaussian kernel filter)
- **sig (float)** (variance of 1d gaussian kernel.)
- **kernel (tuple)** (kernel size to use for smoothing)

**Returns:** result (1d numpy array)

**Return type:** smoothed signal

`moseq2_pca.util.gaussian_kernel1d` (n=None, sig=3)  
Get 1D gaussian kernel.

**Parameters:**

- **n (int)** (number of points to use.)
- **sig (int)** (variance of kernel to use.)

**Returns:** kernel (1d array)

**Return type:** 1D numpy kernel.

`moseq2_pca.util.get_changepoints` (scores, k=5, sigma=3, peak\_height=0.5, peak\_neighbors=1, baseline=True, timestamps=None)

Compute changepoints distribution and CP Curve.

**Parameters:**

- **scores (3D numpy array)** (*nframes \* r \* c*)
- **k (int)** (*klags - Lag to use for derivative calculation.*)
- **sigma (int)** (*Standard deviation of gaussian smoothing filter.*)
- **peak\_height (float)** (*user-defined peak Changepoint length.*)
- **peak\_neighbors (int)** (*number of peaks in the CP curve.*)
- **baseline (bool)** (*normalize data.*)
- **timestamps (array)** (*loaded timestamps.*)

**Returns:** **cps (2D numpy array)** (*array of values for CP curve*) **normed\_df (1D numpy array)** (*array of values for bar plot*)

`moseq2_pca.util.get_metadata_path` (h5file)

Return path within h5 file that contains the kinect extraction metadata.

**Parameters:** **h5file (str)** (*path to h5 file.*)

**Returns:** **(str)**

**Return type:** path to acquisition metadata within h5 file.

`moseq2_pca.util.get_rps` (frames, rps=600, normalize=True)

Get random projections of frames.

**Parameters:**

- **frames (2D or 3D numpy array)** (*Frames to get dimensions from.*)
- **rps (int)** (*Number of random projections.*)
- **normalize (bool)** (*indicates whether to normalize frames.*)

**Returns:** **rproj (2D or 3D numpy array)**

**Return type:** Computed random projections with same shape as frames

`moseq2_pca.util.get_timestamp_path` (h5file)

Return path within h5 file that contains the kinect timestamps

**Parameters:** **h5file (str)** (*path to h5 file.*)

**Returns:** **(str)**

**Return type:** path to metadata timestamps within h5 file

`moseq2_pca.util.initialize_dask` (nworkers=50, processes=1, memory='4GB', cores=1, wall\_time='01:00:00', queue='debug', local\_processes=False, cluster\_type='local', scheduler='distributed', timeout=10, cache\_path='/Users/aymanzeine/moseq2\_pca', \*\*kwargs)

Initialize dask client, cluster, workers, etc.

**Parameters:**

- **nworkers (int)** (*number of dask workers to initialize*)
- **processes (int)** (*number of processes per worker*)
- **memory (str)** (*amount of memory to allocate to dask cluster*)
- **cores (int)** (*number of cores to use.*)
- **wall\_time (str)** (*amount of time to allow program to run*)
- **queue (str)** (*logging mode*)
- **local\_processes (bool)** (*indicate whether the processes are local*)
- **cluster\_type (str)** (*indicate what cluster to use*)
- **scheduler (str)** (*indicate what scheduler to use*)
- **timeout (int)** (*number of worker timeouts to allow*)
- **cache\_path (str or Pathlike)** (*path to store cached data*)
- **kwargs** (*extra keyword arguments*)

**Returns:** **client** (**dask Client**) (*initialized Client*) **cluster** (**dask Cluster**) (*initialized Cluster*) **workers** (**dask Workers**) (*initialized workers or None if cluster\_type = 'local'*) **cache** (**dask Chest**) (*initialized Chest (cache) object pointing to given cache path*)

`moseq2_pca.util.insert_nans` (timestamps, data, fps=30)

Fills NaN values with 0 in timestamps.

**Parameters:**

- **timestamps (1D array)** (*timestamp time-strs*)
- **data (1D array)** (*timestamp values*)
- **fps (int)** (*frames per second*)

**Returns:** **filled\_data (1D array)** (*filled missing timestamp values.*) **data\_idx (1D array)** (*indices of inserted 0s*) **filled\_timestamps (1D array)** (*filled timestamp-strs*)

`moseq2_pca.util.read_yaml` (yaml\_file)

Reads yaml file and returns dictionary representation of file contents.

**Parameters:** **yaml\_file (str)** (*path to yaml file*)

**Returns:** **return\_dict (dict)**

**Return type:** dict of yaml file contents

`moseq2_pca.util.recursive_find_h5s` (root\_dir='/Users/aymanzeine/Desktop/moseq/moseq2-pca/docs', ext='.h5', yaml\_string='{}.yaml')

Recursively find h5 files, along with yaml files with the same basename

**Parameters:**

- **root\_dir (str or os.Pathlike)** (*path to directory to start recursive search*)
- **ext (str)** (*extension to search for, e.g. .h5*)
- **yaml\_string (str)** (*a format to use to name yaml files*)

**Returns:** **h5s (list)** (*list of h5 file paths*) **dicts (list)** (*list of dicts containing metadata file contents*) **yamls (list)** (*list of yaml file paths*)

`moseq2_pca.util.recursively_load_dict_contents_from_group` (h5file, path)

Reads all contents from h5 and returns them in a nested dict object.

**Parameters:**

- **h5file (str)** (*path to h5 file*)
- **path (str)** (*path to group within h5 file*)

**Returns:** **ans (dict)**

**Return type:** dictionary of all h5 group contents

`moseq2_pca.util.select_strel` (string='e', size=10, 10)

Selects Structuring Element Shape

**Parameters:**

- **string (str)** (*e for Ellipse, r for Rectangle*)
- **size (tuple)** (*size of StructuringElement*)

**Returns:** **strel (cv2.StructuringElement)**

**Return type:** returned StructuringElement with specified size.

`moseq2_pca.util.shutdown_dask` (scheduler, workers=None)

Graceful

shutdown

dask

scheduler.

source:

<https://github.com/dask/distributed/issues/1703#issuecomment-361291492>

**Parameters:** **scheduler (dask Scheduler)** (*scheduler to shutdown.*)

**Returns:**

**Return type:** None

## Visualization Module

`moseq2_pca.viz.changepoint_dist` (cps, headless=False)

Creates bar plot describing computed Changepoint Distribution.

**Parameters:**

- **cps (np.ndarray)** (*changepoints to graph*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.display_components` (components, cmap='gray', headless=False)

Creates grid of computed Principal Components.

**Parameters:**

- **components (2D np.ndarray)** (*components to graph*)
- **cmap (str)** (*color map to use*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)** (*figure to save/graph*) **ax (plt.ax)** (*figure axis variable*)

`moseq2_pca.viz.scrree_plot` (explained\_variance\_ratio, headless=False)

Creates Scree plot describing principal components.

**Parameters:**

- **explained\_variance\_ratio (1D np.array)** (*explained variance ratio of each principal component*)
- **headless (bool)** (*trim first element in PC list*)

**Returns:** **plt (plt.figure)**

**Return type:** figure to save/graph

## Subpackages

### `moseq2_pca.helpers` package

#### Helpers - Data Module

`moseq2_pca.helpers.data.get_pca_yaml_data` (pca\_yaml)

Reads PCA yaml file and returns metadata

**Parameters:** **pca\_yaml (str)** (*path to pca.yaml*)

**Returns:** **use\_fft (bool)** (*indicates whether to use FFT*) **clean\_params (dict)** (*dict of image filtering parameters*) **mask\_params (dict)** (*dict of mask parameters*) **missing\_data (bool)** (*indicates whether to use mask\_params*)

`moseq2_pca.helpers.data.load_pcs_for_cp` (pca\_file\_components, config\_data)

Load computed Principal Components for Model-free Changepoint Analysis.

**Parameters:**

- **pca\_file\_components (str)** (*path to pca h5 file to read PCs*)
- **config\_data (dict)** (*config parameters*)

**Returns:** **pca\_components (str)** (*path to pca components*) **changepoint\_params (dict)** (*dict of relevant changepoint parameters*) **cluster (dask Cluster)** (*Dask Cluster object.*) **client (dask Client)** (*Dask Client Object*) **workers (dask Workers)** (*initialized workers or None if cluster\_type = 'local'*) **cache (dask Chest)** (*initialized Chest (cache) object pointing to given cache path*) **missing\_data (bool)** (*Indicates whether to use mask\_params*) **mask\_params (dict)** (*Mask parameters to use when computing CPs*)

`moseq2_pca.helpers.data.setup_cp_command` (input\_dir, config\_data, output\_dir, output\_file, output\_directory=None)

Helper function for changepoints\_wrapper to perform data-path existence checks.

**Parameters:**

- **input\_dir (int)** (path to directory containing all h5+yaml files)
- **config\_data (dict)** (dict of relevant PCA parameters (image filtering etc.))
- **output\_dir (str)** (path to directory to store PCA data)
- **output\_file (str)** (pca model filename)
- **output\_directory (str)** (alternative output\_dir)

**Returns:** **config\_data (dict)** (updated config\_data dict with the proper paths) **pca\_file\_components (str)** (path to trained pca file) **pca\_file\_scores (str)** (path to pca\_scores file) **h5s (list)** (list of relevant pca h5 files) **yamls (list)** (list of relevant pca metadata yaml files) **save\_file (str)** (path to save changepoints)

## Helpers - Wrapper Module

`moseq2_pca.helpers.wrappers.apply_pca_wrapper` (input\_dir, config\_data, output\_dir, output\_file, output\_directory=None, gui=False)

Wrapper function to obtain PCA Scores.

**Parameters:**

- **input\_dir (int)** (path to directory containing all h5+yaml files)
- **config\_data (dict)** (dict of relevant PCA parameters (image filtering etc.))
- **output\_dir (str)** (path to directory to store PCA data)
- **output\_file (str)** (pca model filename)
- **output\_directory (str)** (alternative output\_dir)
- **gui (bool)** (indicate GUI is running)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.compute_changepoints_wrapper` (input\_dir, config\_data, output\_dir, output\_file, gui=False, output\_directory=None)

Wrapper function to compute model-free (PCA based) Changepoints.

**Parameters:**

- **input\_dir (int)** (path to directory containing all h5+yaml files)
- **config\_data (dict)** (dict of relevant PCA parameters (image filtering etc.))
- **output\_dir (str)** (path to directory to store PCA data)
- **output\_file (str)** (pca model filename)
- **output\_directory (str)** (alternative output\_dir)
- **gui (bool)** (indicate GUI is running)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

`moseq2_pca.helpers.wrappers.train_pca_wrapper` (input\_dir, config\_data, output\_dir, output\_file, output\_directory=None, gui=False)

Wrapper function to train PCA.

**Parameters:**

- **input\_dir (int)** (path to directory containing all h5+yaml files)
- **config\_data (dict)** (dict of relevant PCA parameters (image filtering etc.))
- **output\_dir (str)** (path to directory to store PCA data)
- **output\_file (str)** (pca model filename)
- **output\_directory (str)** (alternative output\_dir)
- **gui (bool)** (indicate GUI is running)

**Returns:** **config\_data (dict)**

**Return type:** updated config\_data variable to write back in GUI API

## **moseq2\_pca.pca package**

### **PCA - Utilities Module**

`moseq2_pca.pca.util.apply_pca_dask` (pca\_components, h5s, yamls, use\_fft, clean\_params, save\_file, chunk\_size, mask\_params, missing\_data, client, fps=30, h5\_path='/frames', h5\_mask\_path='/frames\_mask')  
"Apply" trained PCA on input frame data to obtain PCA Scores using Distributed Dask cluster.

**Parameters:**

- **pca\_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use\_fft (bool)** (indicate whether to use 2D-FFT)
- **clean\_params (dict)** (dictionary containing filtering options)
- **save\_file (str)** (path to pca\_scores filename to save)
- **chunk\_size (int)** (size of chunks to process)
- **mask\_params (dict)** (dictionary of masking parameters (if missing data))
- **missing\_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)
- **h5\_path (str)** (path to frames within selected h5 file (default: '/frames'))
- **h5\_mask\_path (str)** (path to masked frames within selected h5 file (default: '/frames\_mask'))

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.apply_pca_local` (pca\_components, h5s, yamls, use\_fft, clean\_params, save\_file, chunk\_size, mask\_params, missing\_data, fps=30, h5\_path='/frames', h5\_mask\_path='/frames\_mask')  
"Apply" trained PCA on input frame data to obtain PCA Scores using local cluster/platform.

**Parameters:**

- **pca\_components (np.array)** (array of computed Principal Components)
- **h5s (list)** (list of h5 files)
- **yamls (list)** (list of yaml files)
- **use\_fft (bool)** (indicate whether to use 2D-FFT)
- **clean\_params (dict)** (dictionary containing filtering options)
- **save\_file (str)** (path to pca\_scores filename to save)
- **chunk\_size (int)** (size of chunks to process)
- **mask\_params (dict)** (dictionary of masking parameters (if missing data))
- **missing\_data (bool)** (indicates whether to use mask arrays.)
- **fps (int)** (frames per second)
- **h5\_path (str)** (path to frames within selected h5 file (default: '/frames'))
- **h5\_mask\_path (str)** (path to masked frames within selected h5 file (default: '/frames\_mask'))

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.compute_explained_variance` (s, nsamples, total\_var)

Computes the explained variance and explained variance ratio contributed by each computed Principal Component.

**Parameters:**

- **s (1d array)** (*computed singular values.*)
- **nsamples (int)** (*number of included samples.*)
- **total\_var (float)** (*total variance captured by principal components.*)

**Returns:** **explained\_variance (1d-array)** (*list of floats denoting the explained variance per PC.*)  
**explained\_variance\_ratio (1d-array)** (*list of floats denoting the explained variance ratios per PC.*)

`moseq2_pca.pca.util.compute_svd(dask_array, mean, rank, iters, missing_data, mask, recon_pcs, min_height, max_height, client, gui=False)`

Runs Singular Vector Decomposition on the inputted frames of shape (nframes, nfeatures). Data is centered by subtracting it by the mean value of the data. If `missing_data == True`, It will iteratively recompute the svd on the mean-centered data to reconstruct the PCs from the missing data until it converges.

**Parameters:**

- **dask\_array (dask 2d-array)** (*Reshaped input data array of shape (nframes x nfeatures)*)
- **mean (1d array)** (*Means of each row in dask\_array.*)
- **rank (int)** (*Rank of the desired thin SVD decomposition.*)
- **iters (int)** (*Number of SVD iterations*)
- **missing\_data (bool)** (*Indicates whether to compute SVD with a masked array*)
- **mask (dask 2d-array)** (*None if missing\_data == False, else mask array of shape dask\_array*)
- **recon\_pcs (int)** (*Number of PCs to reconstruct for missing data.*)
- **min\_height (int)** (*Minimum height of mouse above the ground, used to filter reconstructed PCs.*)
- **max\_height (int)** (*Maximum height of mouse above the ground, used to filter reconstructed PCs.*)
- **client (dask Client)** (*Dask client to process batches.*)
- **gui (bool)** (*Indicates to dask to show a progress bar in Jupyter*)

**Returns:** **s (1d array)** (*computed singular values (eigen-values).*) **v (2d array)** (*computed principal components (eigen-vectors).*) **mean (1d array)** (*updated mean of dask array if missing\_data == True.*) **total\_var (float)** (*total variance captured by principal components.*)

`moseq2_pca.pca.util.copy_metadata_to_scores(f, f_scores, uuid)`

Copies metadata from individual session extract h5 files to the PCA scores h5 file.

**Parameters:**

- **f (read-open h5py File)** (*open "results\_00.h5" h5py.File object in read-mode*)
- **f\_scores (read-open h5py File)** (*open "pca\_scores.h5" h5py.File object in read-mode*)
- **uuid (str)** (*uuid of inputted session h5 "f".*)

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.get_changepoints_dask(changepoint_params, pca_components, h5s, yamls, save_file, chunk_size, mask_params, missing_data, client, fps=30, pca_scores=None, progress_bar=False, h5_path='/frames', h5_mask_path='/frames_mask')`

Computes model-free changepoints using PCs and PC Scores on distributed dask cluster.



**Parameters:**

- **changepoint\_params (dict)** (*dict of changepoint parameters*)
- **pca\_components (np.array)** (*computed principal components*)
- **h5s (list)** (*list of h5 files*)
- **yamls (list)** (*list of yaml files*)
- **save\_file (str)** (*path to save changepoint files*)
- **chunk\_size (int)** (*size of chunks to process in dask.*)
- **mask\_params (dict)** (*dict of missing\_data mask parameters.*)
- **missing\_data (bool)** (*indicate whether to use mask\_params*)
- **client (dask Client)** (*initialized Dask Client object*)
- **fps (int)** (*frames per second*)
- **pca\_scores (np.array)** (*computed principal component scores*)
- **progress\_bar (bool)** (*display progress bar*)
- **h5\_path (str)** (*path to frames within selected h5 file (default: '/frames')*)
- **h5\_mask\_path (str)** (*path to masked frames within selected h5 file (default: '/frames\_mask')*)

**Returns:**

**Return type:** None

`moseq2_pca.pca.util.get_timestamps` (f, frames, fps=30)

Reads the timestamps from a given h5 file.

**Parameters:**

- **f (read-open h5py File)** (*open "results\_00.h5" h5py.File object in read-mode*)
- **frames (3d-array)** (*list of 2d frames contained in opened h5 File.*)
- **fps (int)** (*frames per second.*)

**Returns:** timestamps (1d array)

**Return type:** array of timestamps for inputted frames variable

`moseq2_pca.pca.util.mask_data` (original\_data, mask, new\_data)

Create a mask subregion given a boolean mask if missing data flag is used.

**Parameters:**

- **original\_data (3d np.ndarray)** (*input frames*)
- **mask (3d boolean np.ndarray)** (*mask array*)
- **new\_data (3d np.ndarray)** (*frames to use*)

**Returns:** output (3d np.ndarray)

**Return type:** masked data array

`moseq2_pca.pca.util.train_pca_dask` (dask\_array, clean\_params, use\_fft, rank, cluster\_type, client, workers, cache, mask=None, iters=10, recon\_pcs=10, min\_height=10, max\_height=100, gui=False)

Train PCA using dask arrays.

**Parameters:**

- **dask\_array (dask array)** (*chunked frames to train PCA*)
- **clean\_params (dict)** (*dictionary containing filtering parameters*)
- **use\_fft (bool)** (*indicates whether to use 2d-FFT on images.*)
- **rank (int)** (*Matrix rank to use*)
- **cluster\_type (str)** (*indicates which cluster to use.*)
- **client (Dask.Client)** (*client object to execute dask operations*)
- **workers (int)** (*number of dask workers*)
- **cache (str)** (*path to cache directory*)
- **mask (dask array)** (*dask array of masked data if missing\_data parameter==True*)
- **iters (int)** (*number of SVD iterations*)
- **recon\_pcs (int)** (*number of PCs to reconstruct. (if missing\_data = True)*)
- **min\_height (int)** (*minimum mouse height from floor in (mm)*)
- **max\_height (int)** (*maximum mouse height from floor in (mm)*)
- **gui (bool)** (*Indicates to dask to show a progress bar in Jupyter*)

**Returns:** **output\_dict (dict)**

**Return type:** dictionary containing PCA training results.

## Indices and tables

- **genindex**
- **modindex**
- **search**

# Index

## Symbols

	<code>--h5-mask-path &lt;h5_mask_path&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option
		<code>moseq2-pca-compute-changepoints</code>	command line option
<code>--chunk-size &lt;chunk_size&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--cluster-type &lt;cluster_type&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--config-file &lt;config_file&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--cores &lt;cores&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
	<code>--mask-height-threshold &lt;mask_height_threshold&gt;</code>	<code>moseq2-pca-train-pca</code>	command line option
	<code>--mask-threshold &lt;mask_threshold&gt;</code>	<code>moseq2-pca-train-pca</code>	command line option
	<code>--max-height &lt;max_height&gt;</code>	<code>moseq2-pca-train-pca</code>	command line option
<code>--dask-cache-path &lt;dask_cache_path&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--detrend-window &lt;detrend_window&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
<code>--dims &lt;dims&gt;</code>	<code>moseq2-pca-compute-changepoints</code>	command line option	
<code>--fill-gaps &lt;fill_gaps&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
<code>--fps &lt;fps&gt;</code>	<code>moseq2-pca-apply-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option
	<code>moseq2-pca-compute-changepoints</code>	command line option	
	<code>moseq2-pca-clip-scores</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--gaussfilter-space &lt;gaussfilter_space&gt;</code>	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-compute-changepoints</code> command line option
<code>--gaussfilter-time &lt;gaussfilter_time&gt;</code>	<code>moseq2-pca-train-pca</code>	command line option	<code>moseq2-pca-train-pca</code> command line option

--output-dir <output_dir>	moseq2-pca-apply-pca line option	command	-use-fft	moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option		-visualize-results <visualize_results>	moseq2-pca-compute-char command line option
	moseq2-pca-train-pca line option	command		moseq2-pca-train-pca line option
--output-file <output_file>	moseq2-pca-apply-pca line option	command	--wall-time <wall_time>	moseq2-pca-apply-pca line option
	moseq2-pca-compute-changepoints command line option			moseq2-pca-compute-changepoints command line option
	moseq2-pca-train-pca line option	command		moseq2-pca-train-pca command line option
--pca-file <pca_file>	moseq2-pca-apply-pca command line option	-c		moseq2-pca-apply-pca command line option
-e-components <pca_file_components>	moseq2-pca-compute-changepoints command line option			moseq2-pca-compute-changepoints command line option
-a-file-scores <pca_file_scores>	moseq2-pca-compute-changepoints command line option			moseq2-pca-train-pca command line option
--pca-path <pca_path>	moseq2-pca-apply-pca line option	command	-d	moseq2-pca-apply-pca command line option
	moseq2-pca-compute-changepoints command line option			moseq2-pca-compute-changepoints command line option
--processes <processes>	moseq2-pca-apply-pca line option	command		moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option	-i		moseq2-pca-apply-pca command line option
	moseq2-pca-train-pca line option	command		moseq2-pca-compute-changepoints command line option
--queue <queue>	moseq2-pca-apply-pca line option	command		moseq2-pca-train-pca command line option
	moseq2-pca-compute-changepoints command line option	-k		moseq2-pca-compute-changepoints command line option
	moseq2-pca-train-pca line option	command	-m	moseq2-pca-apply-pca command line option
--rank <rank>	moseq2-pca-train-pca command line option			moseq2-pca-compute-changepoints command line option
--recon-pcs <recon_pcs>	moseq2-pca-train-pca command line option			moseq2-pca-train-pca command line option
--sigma <sigma>	moseq2-pca-compute-changepoints command line option	-n		moseq2-pca-apply-pca command line option
--tailfilter-shape <tailfilter_shape>	moseq2-pca-train-pca command line option			moseq2-pca-compute-changepoints command line option
--tailfilter-size <tailfilter_size>	moseq2-pca-train-pca command line option			moseq2-pca-train-pca command line option
--threshold <threshold>	moseq2-pca-compute-changepoints command line option	-o		moseq2-pca-apply-pca command line option
--timeout <timeout>	moseq2-pca-apply-pca line option	command		moseq2-pca-compute-changepoints command line option
	moseq2-pca-compute-changepoints command line option			moseq2-pca-train-pca command line option
	moseq2-pca-train-pca line option	command		

`-p`      `moseq2-pca-apply-pca`    command line option

`moseq2-pca-compute-changepoints` command line option

`moseq2-pca-train-pca`    command line option

`-q`      `moseq2-pca-apply-pca`    command line option

`moseq2-pca-compute-changepoints` command line option

`moseq2-pca-train-pca`    command line option

`-s`      `moseq2-pca-compute-changepoints` command line option

`-w`      `moseq2-pca-apply-pca`    command line option

`moseq2-pca-compute-changepoints` command line option

`moseq2-pca-train-pca`    command line option

## A

`apply_pca_command()` (in module `moseq2_pca.gui`)

`apply_pca_dask()` (in module `moseq2_pca.pca.util`)

`apply_pca_local()` (in module `moseq2_pca.pca.util`)

`apply_pca_wrapper()` (in module `moseq2_pca.helpers.wrappers`)

## C

`changepoint_dist()` (in module `moseq2_pca.viz`)

`clean_frames()` (in module `moseq2_pca.util`)

### CLIP\_SAMPLES

`moseq2-pca-clip-scores` command line option

`command_with_config()` (in module `moseq2_pca.util`)

`compute_changepoints_command()` (in module `moseq2_pca.gui`)

`compute_changepoints_wrapper()` (in module `moseq2_pca.helpers.wrappers`)

`compute_explained_variance()` (in module `moseq2_pca.pca.util`)

`compute_svd()` (in module `moseq2_pca.pca.util`)

`copy_metadatas_to_scores()` (in module `moseq2_pca.pca.util`)

## D

`display_components()` (in module `moseq2_pca.viz`)

## G

`gauss_smooth()` (in module `moseq2_pca.util`)

`gaussian_kernel1d()` (in module `moseq2_pca.util`)

`get_changepoints()` (in module `moseq2_pca.util`)

`get_changepoints_dask()` (in module `moseq2_pca.pca.util`)

`get_metadata_path()` (in module `moseq2_pca.util`)

`get_pca_yaml_data()` (in module `moseq2_pca.helpers.data`)

`get_rps()` (in module `moseq2_pca.util`)

`get_timestamp_path()` (in module `moseq2_pca.util`)

`get_timestamps()` (in module `moseq2_pca.pca.util`)

## I

`initialize_dask()` (in module `moseq2_pca.util`)

`insert_nans()` (in module `moseq2_pca.util`)

## L

`load_pcs_for_cp()` (in module `moseq2_pca.helpers.data`)

## M

`mask_data()` (in module `moseq2_pca.pca.util`)

### module

`moseq2_pca.gui`

`moseq2_pca.helpers.data`

`moseq2_pca.helpers.wrappers`

`moseq2_pca.pca.util`

`moseq2_pca.util`

`moseq2_pca.viz`

### moseq2-pca-apply-pca command line option

`--chunk-size <chunk_size>`

`--cluster-type <cluster_type>`

`--config-file <config_file>`

`--cores <cores>`

`--dask-cache-path <dask_cache_path>`

`--detrend-window <detrend_window>`

`--fill-gaps <fill_gaps>`

`--fps <fps>`

`--h5-mask-path <h5_mask_path>`

`--h5-path <h5_path>`

`--input-dir <input_dir>`

`--memory <memory>`

`--nworkers <nworkers>`

`--output-dir <output_dir>`

`--output-file <output_file>`

--pca-file <pca\_file>  
--pca-path <pca\_path>  
--processes <processes>  
--queue <queue>  
--timeout <timeout>  
--wall-time <wall\_time>  
-c  
-d  
-i  
-m  
-n  
-o  
-p  
-q  
-w

#### **moseq2-pca-clip-scores command line option**

--from-end  
CLIP\_SAMPLES  
PCA\_FILE

#### **moseq2-pca-compute-changepoints command line option**

--chunk-size <chunk\_size>  
--cluster-type <cluster\_type>  
--config-file <config\_file>  
--cores <cores>  
--dask-cache-path <dask\_cache\_path>  
--dims <dims>  
--fps <fps>  
--h5-mask-path <h5\_mask\_path>  
--h5-path <h5\_path>  
--input-dir <input\_dir>  
--klags <klags>  
--memory <memory>  
--neighbors <neighbors>  
--nworkers <nworkers>  
--output-dir <output\_dir>  
--output-file <output\_file>  
--pca-file-components <pca\_file\_components>  
--pca-file-scores <pca\_file\_scores>  
--pca-path <pca\_path>  
--processes <processes>  
--queue <queue>  
--sigma <sigma>

--threshold <threshold>  
--timeout <timeout>  
--visualize-results <visualize\_results>  
--wall-time <wall\_time>  
-c  
-d  
-i  
-k  
-m  
-n  
-o  
-p  
-q  
-s  
-w

#### **moseq2-pca-train-pca command line option**

--chunk-size <chunk\_size>  
--cluster-type <cluster\_type>  
--config-file <config\_file>  
--cores <cores>  
--dask-cache-path <dask\_cache\_path>  
--gaussfilter-space <gaussfilter\_space>  
--gaussfilter-time <gaussfilter\_time>  
--h5-mask-path <h5\_mask\_path>  
--h5-path <h5\_path>  
--input-dir <input\_dir>  
--local-processes <local\_processes>  
--mask-height-threshold <mask\_height\_threshold>  
--mask-threshold <mask\_threshold>  
--max-height <max\_height>  
--medfilter-space <medfilter\_space>  
--medfilter-time <medfilter\_time>  
--memory <memory>  
--min-height <min\_height>  
--missing-data  
--missing-data-iters <missing\_data\_iters>  
--nworkers <nworkers>  
--output-dir <output\_dir>  
--output-file <output\_file>  
--processes <processes>  
--queue <queue>  
--rank <rank>

```

--recon-pcs <recon_pcs>
--tailfilter-shape <tailfilter_shape>
--tailfilter-size <tailfilter_size>
--timeout <timeout>
--use-fft
--visualize-results <visualize_results>
--wall-time <wall_time>

-c
-d
-i
-m
-n
-o
-p
-q
-w

```

### **moseq2\_pca.gui**

module

### **moseq2\_pca.helpers.data**

module

### **moseq2\_pca.helpers.wrappers**

module

### **moseq2\_pca.pca.util**

module

### **moseq2\_pca.util**

module

### **moseq2\_pca.viz**

module

## ***T***

train\_pca\_command() (in module moseq2\_pca.gui)

train\_pca\_dask() (in module moseq2\_pca.pca.util)

train\_pca\_wrapper() (in module  
moseq2\_pca.helpers.wrappers)

## ***P***

### **PCA\_FILE**

moseq2-pca-clip-scores command line option

## ***R***

read\_yaml() (in module moseq2\_pca.util)

recursive\_find\_h5s() (in module moseq2\_pca.util)

recursively\_load\_dict\_contents\_from\_group() (in  
module moseq2\_pca.util)

## ***S***

scree\_plot() (in module moseq2\_pca.viz)

select\_strel() (in module moseq2\_pca.util)

setup\_cp\_command() (in module  
moseq2\_pca.helpers.data)

shutdown\_dask() (in module moseq2\_pca.util)





# Python Module Index

## *m*

[moseq2\\_pca](#)

[moseq2\\_pca.gui](#)

[moseq2\\_pca.helpers.data](#)

[moseq2\\_pca.helpers.wrappers](#)

[moseq2\\_pca.pca.util](#)

[moseq2\\_pca.util](#)

[moseq2\\_pca.viz](#)